I3029P | Ejercitador de piernas
جهاز تمرين الساقين
Beintrainingsgerät
Leg exerciser
Exerciseur des jambes
Pedaliera
Pedaleira de pernas
Тренажер для ног

(ES)

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INDICATIONS:

- For strengthening the leg muscles. It can also be used to exercise the arms.
- This exerciser is intended for home use, not for sport or by professionals. It is suitable for people with a low level of physical activity.

INSTRUCTIONS FOR USE:

- Before embarking on any exercise programme, you should consult with your doctor. Remember to breathe properly, exercise at your own pace, and don't exert yourself to an extent that could harm you.
- Read and follow these instructions carefully before using the exerciser. If you have any questions, always consult your doctor or the supplier you purchased it from before use.
- Always follow these general instructions for use. If you notice any adverse effects, stop exercising and contact your doctor.
- Before each use, carefully check the condition of all the components of the exerciser. If any component is not in a serviceable condition, it should be replaced before use. Also check, before each use, that it can rotate freely and without obstacles, avoiding anything that may limit the mobility of the exerciser's components.
- With the leg exerciser on the floor, you can pedal while sat down to strengthen your leg muscles. By placing it on a table, you can also exercise your arm muscles.
- The level of effort is regulated by the adjustment device (resistance screw) located in the centre of the exerciser. Turn it clockwise to increase resistance and counter-clockwise to decrease resistance.
- It can be easily stored under a bed or a chair. The structure is made of double-coated steel, and the feet are made from non-slip rubber, as are the pedals and their straps. The exerciser is lightweight and requires no assembly. Clean it with a damp cloth. Do not use cleaners with solvents.
- Position the rubber feet so that the exerciser rests safely and stably on the chosen surface.
- The resistance screw is attached to one of the pedals. Release it and insert it into the exerciser.

WARNINGS:

- This exerciser must not be used for any other purpose than that prescribed. Do not let children play or climb on the exerciser.
- **Do not stand on the pedals. Do not lean on the exerciser.** Always use it while seated.
- To avoid the risk of fingers or hands being trapped, do not handle the internal parts of the exerciser while in use.
- Make sure there's at least one metre of free space in all directions when using the exerciser.
- Make sure the exerciser is resting on a flat, horizontal, and solid surface. Take the necessary precautions so that it does not slip during use.
- Wear appropriate clothing and footwear when you exercise. It should not be used with bare feet or loose footwear (flip flops, etc.). Do not wear overly loose clothing that could get caught during exercise.
- Due to the friction caused by the resistance applied to the exerciser, some of the metal parts may be hot after use. Pay attention when handling it.

USEFUL LIFE:

If the conditions of use are respected, its theoretical useful life is one year.

NOTE: Any serious incident related to the product must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient resides.

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PRIM, S.A
Homecare Enterprise Co., Ltd