MANUAL DE INSTRUCCIONES

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LIBRA

Manual de instrucciones **ES**

Instruction manual **EN**

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Manuale D'Istruzioni **IT**

Anweisungen **DE**

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Bruksanvisning **SWE**

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NOTICE

Notice to the user and/or patient: any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is established.

1. INDICATIONS

1. The device should only be used by one user.
2. A slower speed should always be used when performing the following actions: ascending, descending or traversing a slope or inclined terrain, as well as on uneven terrain, ramps and soft surfaces.
3. Keep away from flammable objects while the device is charging.

2. CONTRAINDICATIONS

1. Never attempt to ascend or descend a slope on slippery surfaces or where there is a danger of skidding.
2. Do not drive on a slope steeper than 6 degrees.
3. Avoid shifting the centre of gravity when turning left or right.
4. Direct exposure to rain or moisture will cause the scooter to malfunction and may cause the scooter to rust.
5. Never put the scooter in neutral when on slopes.
6. Do not attempt to modify or tamper with the mechanism or structure of the scooter.
7. Do not use the scooter on stairs or escalators.
8. Do not use the device while it is charging.

3. IMPLEMENTATION

You can make the following adjustments to the scooter to increase your comfort while riding:

- Seat height and distance to the handlebars.
- Steering column.
- Legroom.

1. To start using the scooter, carry out the following steps:

- Remove the rear casing and connect the batteries.
- Check that the immobiliser lever is in D mode.
- Check that the speed control is in low speed (turtle) mode.

2. Stand on the scooter and turn the key. The battery indicator will light up. Confirm that the self-test indicator is not flashing.

3. By placing your hands on the handlebars, you should be able to comfortably reach the control levers, which are operated with your thumbs. The right-hand lever activates the scooter to move forwards and the left-hand lever moves the scooter backwards. If you release both levers, the scooter will stop.

4. To turn the scooter, you must turn the steering column in the direction you want to go.

5. Start practising in obstacle-free areas. Start the scooter at minimum speed and drive backwards and forwards. Try a few turns. As you become more comfortable, you can increase your speed. To do this, turn the control towards the hare icon.
6. When only two of the bars on the battery level LED indicator are lit, you should recharge the battery as soon as possible.
7. If the scooter stops, locate the circuit breaker on the top of the battery. Press it and test the scooter again.
8. When you have finished using the scooter, turn the key to disconnect the scooter before getting off.
9. If you will not be using the scooter until the next day, recharge the battery.

Additional information on operation


When going uphill, you may have to use a higher speed. Remember to slow down again before going downhill.

When going downhill, moderate your speed. To do this, adjust the speed control towards the turtle icon. This allows the user to take more control over the direction of the scooter, as the closer the control is placed to the turtle icon, the slower the scooter will move. However, this model of scooter will not self-accelerate when travelling downhill thanks to the automatic braking system that will be activated if the user tries to travel at too high a speed.

If the scooter stops working and must be moved, move the starter lever to the N position, move the scooter carefully to a safe place and return the lever to the D position.

4. IDENTIFICATION OF THE PARTIES

5. FUNCTIONING OF THE PARTIES

Control panel

- 1. Speed control:** The tortoise icon indicates a slow speed level. The hare icon indicates a fast speed level. Turning this knob will control the total speed that is transferred to the thumb-activated controls.
- 2. Self-diagnostic warning indicator:** A flashing light indicates that there is a problem with the scooter.
- 3. Battery level indicator:** Below this indicator, there are three colour ranges: red, yellow and green. When the light is on and all three colours are displayed, the battery is fully charged. When the light is on, but only shows red, or red areas and yellow areas, the battery needs to be recharged.

Below the control panel

- 1. Throttle:** Pressing the right hand throttle moves the scooter forward. Pressing the left accelerator pedal moves the scooter backwards (this can be reversed if required by your local dealer). Releasing both activates the automatic brake. This is also your throttle. The harder you squeeze them, the faster you go. (Subject to the position of the Rabbit/Turtle control).

Handlebar base

- 1. Adjusting the steering handlebar:** Push the handlebar down to adjust it. Pull up to lock the turn of the steering column at an angle that is comfortable for you.

Under the seat

1. Seat adjustment lever:

- Handlebar clearance adjustment: pull up on the lever on the front of the seat to adjust the handlebar clearance.
- Rotation adjustment: press the lever backwards to rotate the seat; press it forwards to lock it in the desired position.
- Armrest adjustment: Loosen the two thumbscrews to adjust the width between the armrests. Tighten them again to fix them in the desired position.

Rear

1. **Rear cover:** serves to protect the battery from damage.
2. **Anti-tipping wheels:** serve to prevent the scooter from tipping over.
3. **Starter lever:** The starter lever can be used to switch the scooter into manual mode (N) or automatic mode (D). Manual mode allows the user to move the scooter without using the battery, i.e. by pushing it, while automatic mode allows the user to both drive the scooter and keep it stationary when parked. Indicators on the lever will help you determine which mode the scooter is in. Whenever you want to change from manual to automatic mode, you must first switch the scooter off. You must then change mode using the lever. Then switch on the scooter to start using it.

6. MOUNTING THE SCOOTER

1. Raise the steering column using the lever.
2. Tilt the backrest until it is horizontal. Insert the frame tube into the scooter. Make sure that the safety pins have been removed.
3. Fit the two parts of the unit together until the front of the unit stops.
4. Insert the two safety pins into the corresponding holes. Both pins must be fully inserted when the scooter is mounted.
5. Fit the battery, use the Velcro to secure in position and connect the power sockets.
6. Fit the rear cover.
7. Screw on the controller wiring harnesses.
8. Fit the rear compartment cover.
9. Release the seat rotation lever by turning it anticlockwise. Place the seat on the seat post. Retighten the lever by turning it clockwise. Remember to secure the seat control lever before sitting on the scooter.

7. ELECTROMAGNETIC INTERFERENCE WARNINGS

The operation of electric wheelchairs and motorised scooters can be affected by electromagnetic interference (EMI). This is the result of interference from electromagnetic energy (EM) waves emitted from sources such as radio stations, television stations, amateur radio transmissions, two-way radios and mobile phones. Interference (from radio wave sources) could cause the braking system to release unintentionally, cause the scooter to move by itself or to move in an undesirable direction. It could also cause irreparable damage to the scooter's control system.

1. Do not operate portable transceivers (transceivers), such as citizens band (CB) radios. Do not activate personal communication devices, such as mobile phones, while the scooter is on.
2. Be aware of the existence and location of nearby transmitters, such as radio or television stations, and avoid getting too close to them.
3. In the event of unintended movement or brake release, switch the scooter off.
4. Please note that fitting accessories or components to the scooter or modifying the scooter in any way may make it more susceptible to the effect of electromagnetic interference. The immunity level of 20 V/m is generally an affordable and practical level of protection against electromagnetic interference (as of May 1994). Remember that the higher the immunity level, the higher the level of protection. This product has an immunity level of at least 10 V/m.

8. BATTERY CHARGING

The battery should be charged before using the scooter for the first time and should be recharged after each day's use. It is advisable to let the battery charge overnight for the first time and it should be recharged after each day's use. To do this, you will need the scooter and the battery charger.


1. Insert the battery cable into the battery connector.
2. Plug the other end of the cable into a power socket.
3. Switch on the scooter. Normally, the power LED indicator will light up when the power starts to flow.
4. This will start the charging process. During this process, the charging LED will show an orange light. When the light turns green, the battery is fully charged.
5. LED indicator
 - LED (Power) - Green light: Scooter on
 - LED (Charging) - Orange light: Charging activated / Green light: Battery fully charged
6. Troubleshooting loading problems.
 - (A) If the LED (Power) indicator is off:
 - Check that the input voltage (115V/230V) is the same as the voltage that has been set.
 - If the indicator still does not light, check and repair the battery charger.
 - (B) If the LED (Battery) indicator is off:
 - Check that the connection system is working properly.
 - If the battery is fully charged, the LED (Charge) indicator will go out.
 - If the light is still off, the battery may be defective.
 - (C) If the ORANGE light turns GREEN:
 - The battery cannot be charged. You will need to check if it can be repaired.
 - (D) If the ORANGE light turns immediately to GREEN:
 - Check whether the battery is fully charged. If not, it may be defective. You will need to check if it can be repaired.

It is recommended that the batteries are fully charged to avoid damaging the life of the batteries. If the batteries are not fully charged, the scooter may fail.


9. DISMANTLING THE SCOOTER

You do not need any tools to disassemble the scooter. You can disassemble it when you need to transport it.


First, remove the seat. To do this, release the seat rotation lever (A) and lift it up. Carefully lift the scooter's rear compartment cover (B).


Release the wiring harnesses from the controller (C). Remove the scooter rear cover (D).

Remove the Velcro fastening from the battery (E). Disconnect the two battery harnesses (F) and remove the battery (G).


Disconnect the front and rear sections. To do this, hold the seat post (J) with one hand and use the other hand to push the lever towards the rear of the scooter. Lift the seat post to separate the rear frame (H) from the front frame.

10. FUSE REPLACEMENT

1. Remove the screws of the front cover (R), the arrow indicates the direction. (Caution, the screws have different sizes).

2. Then open the fuse box (R1) and replace it with a back-up fuse (R2); note that the fuse must be of the same rating.

11. CLEANING AND MAINTENANCE OF THE SCOOTER

Cleaning the scooter: If the scooter is dirty, use a damp or slightly soapy cloth to clean it. Do not wet it directly or rinse it, as it has electrical parts that should not come into direct contact with water. You can use car polish liquid. All maintenance and repair work on the scooter should be carried out by a specialist. The following areas require regular inspections:

- Tyre pressure should be between 30 and 50 psi.
- Check that all electrical connections are securely connected.

TECHNICAL SPECIFICATIONS	
Overall length	128 cm
Torso width (with armrests)	61 cm
Total height	111 cm
Front wheels	26 cm
Rear wheels	26 cm
Weight	87.9 kg
Maximum speed	8 km/h
Maximum supported weight	135 kg
Ground clearance	4 cm
Degree of inclination	12 degrees
Maximum kerb height	7 cm
Turning radius	152 cm
Brake	Electromechanical
Type of seat	Swivel and sliding seat with head restraints
Seat width	45.5 cm
Engine size	500 W. 3400 r.p.m
Batteries	(2) 12V. 36Ah
Battery pack weight	22.5 kg

Autonomy	30 km
Battery charger	5A external

12. PROBLEM DIAGNOSIS

Count the number of times the blue light on your scooter flashes to find out your fault.

NUMBER OF FLASHES	DESCRIPTION	MEANING
1	Low battery	The batteries are exhausted. - Recharge the battery.
2	Battery depleted	The batteries are exhausted. - Recharge the battery.
3	Battery overload fault	Battery voltage is too high. If the charger is connected, disconnect it. May have been caused by overcharging or overcharging or driving downhill. - When driving downhill, reduce speed. - Check its condition and connections.
4	Engine runtime failure	The engine running time has been exceeded. Switch the scooter off and leave it for a few minutes before switching it on again. - Check the motor, connections and wiring for short circuits. there is a short circuit. - Switch off the controller and allow a few minutes before switching back on.
5	Immobilisation brake failure	Immobiliser brake is activated or faulty: - Make sure that the starter lever is in automatic mode (D). mode (D). - If the flashing light does not come on until the throttle is used, check that there is no short circuit in the brake. - If the light appears as soon as the scooter is switched on, check that there is no open circuit fault. - Check the brake, its connections and wiring.
6	Throttle failure	The throttle is not in neutral when the scooter is switched on. - Make sure the throttle is in neutral mode before starting your scooter. - The throttle may need to be re-adjusted. - Check the throttle connections.

7	Speed failure	Throttle or throttle connections are faulty. - Check wiring for short circuits or open connections. open connections. - Throttle may not have been installed correctly.
8	Engine failure	The motor or its wiring is faulty. - Check the motor, its connections and wiring.
9	Another mistake	The controller may have an internal fault. - Check connections and wiring.

ES - Peso máximo 135 kg
EN - Maximum weight 135 kg
IT - Peso massimo 135 kg
FR - Poids maximum 135 kg
PT - Peso máximo 135 kg
DE - Maximales Gewicht 135 kg
NL - Maximaal gewicht 135 kg
SWE - Högsta vikt 135 kg
PL - Maksymalna waga 135 kg
DK - Maksimal vægt 135 kg

CE **MD** Made in P.R.C

CHIEN TI ENTERPRISE CO., LTD
No. 13, Lane 227, Fu Ying Rd., Hsin Chuang District, New Taipei City,
Taiwan
sales@chienti.com.tw.

MEDNET EC-REP GMBH
Borkstrasse 10, 48163, Muenster, Germany
ecrep@medneteuropa.com

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GRUPO R. QUERALTÓ S.A.