# MANUAL DE INSTRUCCIONES

## ELECTROCARDIÓGRAFO


### ECG1200G

Manual de instrucciones

ES

Instruction manual

EN

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
### ECG1200G

**NOTICE**

**Notice to the user and/or patient; any serious incident that has occurred in relation to the product must be informed to the manufacturer and to the local authority of the Member State in which the user and/or patient is established.**

**INDICATIONS**

- In any type of patient to detect episodes of hypoxaemia difficult to detect clinically.
- Situations of respiratory failure of any aetiology.
- Emergency tracheal intubation and connection to mechanical ventilation.
- As a heart rate monitor in any emergency procedure.

**1. INTRODUCTION**

This product is an electrocardiograph with 12 simultaneous ECG leads that allows recording ECG waveforms with a thermal printing system. Its functions are as follows: recording and displaying ECG waveforms in automatic! manual mode; self-measuring and self-diagnosis of ECG waveform parameters; indication of lead off and paper out; language switching (Chinese/English); AC or DC power supply; rhythm lead can be selected to detect abnormal heart rhythms; database management. High resolution 800\*600 LCD display, 1728 dots wide printing, can print twelve-channel ECG waveforms. operated by keyboard or touch screen.

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- Environmental conditions
  - Usage
    - a) Ambient temperature +5°C~+35°C
    - b) Relative humidity: ≤80%.
    - c) Power supply : AC: 100V~240V, 50/60Hz. Rechargeable lithium battery 14.8V, 3700 mAh
    - d) Atmospheric pressure: 86kPa~106 kPa
  - Transport and storage:
    - a) Ambient temperature: -10°C~55°C
    - b) Relative humidity: ≤95%.
    - c) Atmospheric pressure: 50 kPa~106 kPa
- Input mode: Floating circuit and defibrillation protection
- Lead: 12 standard leads,
- Patient leakage current: <10μA
- Input impedance: ≥50MΩ
- Frequency response: 0.05Hz~150Hz (-3dB)
- Time constant: Time constant >3.2s
- CMRR: >60dB, >100dB (Add filter)
- Interference filter: EMG: 35Hz(-3 dB)
- Recording mode: thermal printing system
- Printing paper specifications: 210 mm (width)\*20 m (height), high speed thermal paper.
- Paper speed:
  - Auto registration: 25 mm/s, 50mm/s, error ±5%.
  - Pace registration: 25 mm/s, 50 mm/s, error ±5%.
  - Manual registration: 5 mm/s, 6.25 mm/s, 10 mm/s, 12.5 mm/s, 25 mm/s, 50 mm/s, 100 mm/s, error ±5%.
- Sensitivity options: 2.5, 5, 10, 20, 40 mm/mV, error ±5%. Standard sensitivity is 10 mm/mV ±0.2 mm/mV.
- Automatic recording: works according to recording format and automatic mode, automatic lead change, automatic measurements and analysis.
- Pace recording: operates according to pacing format and pacing mode, automatic measurements and analysis.
- Manual recording: works according to recording format, manual lead change.
- Measurement parameters: HR, P-R interval, P duration, QRS duration, QRS duration, T duration, Q-T interleave, Q-Tc, QRS axis, T axis, R (VS), S (V1), R (VS)+S(V1)
- Product safety type: Class I, applied parts type CF. With anti-defibrillation circuit.
- Bias voltage: ±300mV.
- Noise level: ≤15 μVp-p.
- Fuse specifications: 2 pcs. AC delayed fuse Φ5\*20 mm
- Dimensions: 340 mm (length)\*320 mm (width)\*85 mm (height).
- Net weight: 3.2 kg.

## 2. SAFETY MEASURES

- The equipment must be connected to the mains with an earthed mains supply.
- If the earthing cable is not integrated, the equipment can only be operated with the internal battery.
- Before changing the fuse, unplug the equipment from the mains.
- The device may only be operated and maintained by qualified personnel.
- Before replacing the fuse, disconnect the equipment from the mains.
- The equipment must be operated and maintained by qualified personnel only.
- Before using the device, the user must read the manual carefully and follow the instructions.
- The equipment is designed in compliance with safety regulations, but under no circumstances should the equipment be left running unattended.
- Before cleaning or disinfecting, the equipment must be switched off and unplugged from the mains.
- Do not use the device in the presence of flammable anaesthetic gases.
- If using the equipment with a cardiac defibrillator or other electrostimulation device, the Ag/AgCl chest electrodes and patient cable must be used with appropriate protection. To prevent the metal electrode from burning the skin, use disposable electrodes whenever the defibrillation duration exceeds 55 seconds. It is recommended that this equipment not be used in conjunction with other electrical stimulation equipment at the same time. If you need to do so, the process should be supervised by a professional technician.
- All accessories connected to the equipment must be Type I and comply with IEC60601-1. Because the total leakage current can harm patients, the leakage current must be monitored.
- Special attention should be paid to the following points when interpreting measurements:
  - P and Q-waves are not always reliable in case of muscle artefacts or AC interference. The same is true for the ST segment and the T wave.
  - Blurred ends of S and T waves may lead to tolerance in the measurement process.
  - If the R wave has been left out due to low voltage of the QRS wave or the drop of a lead, the measured heart rate may deviate greatly from the correct value.
  - Calculation of the axis and identification of the QRS line border are not always reliable in case of low QRS wave voltage.
  - Occasionally, frequent ventricular premature complexes may be identified as the dominant rhythm.
  - Mixed arrhythmias may result in unreliable measurement due to the difficulty in distinguishing the P wave in such a situation.
  - This device is designed to perform interpretation of the ECG tracing immediately after measurement. This interpretation does not contain information about all possible heart problems and sometimes may not agree with the physician's diagnosis. Therefore, the final conclusion about each patient is up to the physician based on the patient's symptoms, ECG results, and other variables.

## 3. FEATURES

- 800\*600 LCD colour screen with touch screen and control keys.
- High resolution thermal printer (8 dots/mm). No adjustment required. Frequency response up to 150 Hz.
- The device allows to record clear and precise waveforms of 12 channels in real time and continuously. Display shows: drift, sensitivity, paper speed, filter status, etc.
- In automatic mode, just press a button to start recording.
- Digital filter design to avoid baseline drift and AC and EMG interference.
- Automatic baseline adjustment, optimisation of recording settings, rhythm leads available.
- Keypad and touch screen control for easy device operation. Status is displayed for your convenience.
- Combination AC/DC power system. This device includes a rechargeable lithium battery, automatic protection circuitry and can be set to automatically shut off if no signal is received.
- This device records 150 ECG waves and prints for 90 minutes continuously.
- The device's memory allows 1000 patient records to be stored, which is useful for data review and statistics.
- Sleek and lightweight design.
- You can record and save data at the touch of a button.
- Digital signal processor of AC, DFT and EMG filters to obtain a high quality electrocardiogram.
- The device has the function of automatic measurement, analysis and diagnosis to facilitate the work of professionals and reduce their workload.
- The device collects signals from 12 ECG channels simultaneously and analyses them with different recording modes such as 12\*1, 6\*2+1, 6\*2, rhythm row 12, row 10, row 8, row 6, etc. There are different formats for reports.
- The device can display 3, 6, or 12 waveforms of leads in real time with different views. There are indications of lead off, lack of paper, etc. In addition, parameters can be automatically detected and reports can be created for the analysis performed. Waveforms can also be recorded in manual or automatic mode and patient data is stored in a database.
- Safety classification: Class I, applied parts type CF defibrillation-proof.
- Degree of protection against liquid ingress: IPX0.
- According to the degree of safety when used in environments where flammable anaesthesia gases are mixed with air: this device is not suitable for use in such environments.
- According to the working mode: continuous operation.

•Symbol explanation:

~AC AC power supply  
 OFF Disconnected/Off  
 ON Connected/On

Equipotential point

Warning

Type CF device, with defibrillation protection function

USB connector

Lead connector

## 4. PARTS DESCRIPTION

### 4.1 Panel elements and views

Paper compartment

Display

Control panel

Front View

Side view

Grounding port


Lead connector

USB Port

Power switch






AC Input

Fuse socket

Bottom view

### 4.2. Buttons description

	Turn the device on/off		Start/Stop sampling
	Select filter		Confirm operation
	Adjust paper speed		Up
	Adjust gain		Down
	Change printing mode		Left
	Start printing		Right
	Backspace		
	Switch input method		
	Number/Uppercase and lowercase/Character: 0/Space/ ,		
	Number/Uppercase and lowercase/Character: 1/ .		
	Number/Uppercase and lowercase/Character: 2/abc/ABC/ :		
	Number/Uppercase and lowercase/Character: 3/def/DEF/;		
	Number/Uppercase and lowercase/Character: 4/ghi/GHI/ *		
	Number/Uppercase and lowercase/Character: 5/jkl/JKL/%		
	Number/Uppercase and lowercase/Character: 6/mno/MNO";		
	Number/Uppercase and lowercase/Character: 7/pqrs/PQRS/#		
	Number/Uppercase and lowercase/Character: 8/tuv/TUV/(		
	Number/Uppercase and lowercase/Character: /wxyz/WXYZ/)		

### 5.3. Light indicators

 The indicator turns green when the equipment is connected to the main electric supply and is switched on. The control light changes colour between green and red when the battery is charging.

 Device is powering on.

## 5. PRECAUTIONS BEFORE USE

- Read this instruction manual carefully before use to ensure safe and effective use of the equipment.
- For the installation and maintenance of the equipment, please note the following:
  - There must not be any high voltage cables, X-ray instruments, ultrasound instruments, electrotherapy instruments etc. near the equipment.
  - Do not use or store the equipment in areas where the air pressure is too high, the temperature and humidity are higher than normal, ventilation is not good, there is too much dust or there are gases containing salt, alkali or chemical agents.
- Place the equipment on a flat surface. When moving the equipment, move it gently. Avoid excessive vibration and shock.
- The AC frequency and voltage must meet the requirements and have sufficient current capacity.
- Please install the device in a place where it is easy to ground it. Avoid contact with patients and patient cables with other conductive materials such as the floor or beds.

## 6. PREPARATION BEFORE USE


- Check that the grounding is correct and that the cable connections are secure.
- Ensure that all electrodes are correctly connected to the patient.
- Ensure that a battery with the correct voltage rating is used.
- Apply the electrode gel one electrode at a time. Remember that two electrodes should never touch to avoid short circuits.
- The power cord should not touch the ECG cable.

## 7. PRECAUTIONS DURING USE

- Continuously monitor the condition of the patient and the equipment.
- The patient and the device may only be connected with the ECG cable.
- Ensure that the patient and equipment do not move during use of the device.
- Switch off the device after use.
- Turn off the equipment and remove the ECG channels gently without pulling them.
- Store the device and accessories properly for the next use.
- Loading of printing paper.
  - The dimensions of the high-speed thermal recording paper are: 210 mm (width)\*20 m (length).
  - Open the paper compartment cover, pull out the paper spindle, and insert the roll paper. Then place them back into the paper compartment.
  - Close the compartment cover. Leave 2 cm of paper outside.

## 8. USE OF THE RECORDING PAPER

- During registration, the device will stop printing in case of paper shortage. The following indication will appear on the screen when paper is missing.

- In order to obtain high quality printouts, it is recommended to use high-sensitivity thermal paper. Other paper may produce smudges and may damage the printing mechanism.
- Protect the thermal paper from heat sources, high temperatures, humidity, exposure to sunlight. It should be stored in a dark, dry and cool place.
- Substances that can cause stains on the paper: gel, glue, organic solvent, etc.
- Substances that may cause the print to fade: Soft PVC folders, plastic, magnetic tapes with plasticiser, fluorescent light, stamping ink, etc.


## 9. PLACEMENT OF THE ELECTRODES

It is recommended to place the electrodes on the thorax first and then on the extremities.

### 9.1 Chest electrodes

Observe the following picture and place the chest electrodes as shown below:

- V1: Fourth intercostal space to the right of the sternum.
- V2: Fourth intercostal space to the left of the sternum.
- V3: Fifth rib between V2 and V4.
- V4: Fifth intercostal space on the left mid-clavicular line.
- V5: Left anterior axillary line, in the horizontal plane of V4.
- V6: Left mid axillary line, in the horizontal plane of V4.


Clean the skin with alcohol in the area where the electrodes are to be placed on the chest. Then apply the ECG cream around 25mm in diameter and up to the edge of the electrodes on the chest and connect the electrodes in positions V1-V6.

**Note:** Please note that it is not allowed for the electrodes to come into contact with each other, as well as the overlapping of the cream from one to the other.

### 9.2. Limb electrodes

The electrodes should be placed on the skin of the hands and feet. Clean all limb electrodes and the areas to which they are to be attached with alcohol before applying the ECG cream. Then attach the electrodes.

**Note:** After inserting the patient cable into the patient cable socket, tighten it firmly.

### 9.3. Electrode Location

Electrode Location	Electrode Code	Socket Number
Right arm	RA/R	9
Left arm	LA/L	10
Left leg	LL/F	11
Right leg	RL/N	14
Chest 1	V1/C1	12
Chest 2	V2/C2	1
Thorax 3	v3/c3	2
Thorax 4	v4/c4	3
Thorax 5	v5/c5	4
Thorax 6	v6/c6	5

## 10. GROUNDING AND MAINS CONNECTION

Make sure that the equipment is switched off and check that it is correctly connected to the earthed socket with a 3-pin connector. The power cable must be correctly inserted into the socket on the equipment. Proper grounding will ensure proper operation and performance of the equipment and prevent AC and electromagnetic wave interference.

## 11. USE OF THE BATTERY

- This equipment has been designed to operate with a rechargeable lithium battery. The display will show the battery status of the battery of 5 bars. The battery will automatically recharge when the equipment is connected to the AC power supply. The total time battery charging time is 4 hours.
- If the battery is charged, the equipment can print continuously for 90 minutes and record 4 hours of data (without printing). - When the battery is discharged, the equipment will automatically turn off.
- The battery should be recharged after use and after 3 months of storage. Otherwise the battery life will be shorter.
- Description of symbols related to the power supply.

N°	Symbol	Description of status
a		Status unknown
b		External power supply
c		Internal power supply
d		Internal power supply, status 3/4
e		Internal power supply, status 1/2
f		Internal power supply, status: 1/4
g		Internal power supply is about to discharge. Must connect AC

Note: When charging, the battery icon changes from (f) to (c).


- When the battery does not recharge or operates less than 10 minutes after being fully charged, the battery must be replaced with a new one.

### WARNING

- Do not connect the "+" and "-" poles of the battery directly with the cable, it may cause fire hazard.
- Protect the battery from heat sources. Risk of explosion.
- Never open or disassemble the battery.

## 12. CONTROL PANEL AND BUTTONS

### 12.1. Main display

Battery level: See 12.4

Buttons:

Input to the sampling interface. Starts automatically, after switching on the equipment.

Entry to the file management interface (modification or deletion of files).

The image of a correct electrode placement.

Time and date settings.

System settings.

Sampling configuration.

Configuration of analysis parameters for each parameter.

Print mode settings.

Manufacturer information.

Quick change: you can quickly change function with the     buttons on the keypad and accept the  button.


## 12.2 Sampling interface

Select the button  in the main menu or press the MENU button  to access the sampling interface.

**Note: Due to settings related to the timing of patient data entry, patient data must be inserted prior to sampling (see 12.3 Insert Patient Data).**

The sampling screen offers several views of the leads: 3 leads, 6 leads, 6\*2 leads, 6\*2 and 12 leads. In the next screen the view is 6\*2 leads.

• Stop Sampling: during sampling, you can press  on the main panel to stop it and return to the main screen.

• Change lead: You can change the current lead waveform and it will be printed in the manual recording.

• Change Lead View: You can change the view between 3-lead, 6-lead, 6\*2-lead and 12-lead.


• Lead-off: In Demo mode, it displays "DEMO ECG", WHILE IN SAMPLING MODE, it displays lead-off information detected.

• Change print mode: by pressing  the display "DEMO ECG" appears. In sampling mode, the indication "lead off" appears.

- Adjust gain (sensitivity): with the button  on the control panel, you can switch between 2.5mm/mV, 5mm/mV, 10mm/mV, 20mm/mV, 40mm/mV.

- Adjust paper speed: With the button  on the control panel, you can adjust the paper speed between 5mm/s, 6.25 mm/s, 10mm/s, 12.5 mm/s, 25 mm/s, 50 mm/s.

• Change filter: with the button  on the control panel, you can change the filter between "no filter", AC, EMG, DFT, AC+EMG, AC+DFT, EMG+DFT, AC+EMG+DFT.

- AC: AC filter.

- EMG: EMG filter.

- DFT: Baseline filter.

• Print/stop printing: With the button  on the main panel, you can start or stop printing.

• Auto mode: After starting the printout, the system will automatically print and save the twelve-lead waveform at the same time. The length will be as set in the print options. Depending on the settings, data and test results will be printed. Printing stops automatically.

• Manual mode: After starting the printout, the user can print different waveforms by changing the lead. This means that the ECG waveform printout in manual mode is not synchronous and cannot be saved. The user must stop the printout by pressing the


button again.

Statuses displayed during the printing process:

Message on the screen	Description
Processing...	In process of printing.
Wait...	In process of finishing printing.
No paper	Thermal paper is missing. After putting paper, you must restart printing.
Print timeout	The connection between the system and the printing sub-system is broken.
ECG timeout	The connection between the system and the sampling sub-system is broken.
Low Power	Battery level is not sufficient to start.


### 12.3 Entering patient data

Depending on the configuration (see 13.8), the user can enter the patient data before or after sampling, it is also possible to leave the data blank. This screen looks like the following one:


Select any field. When pressing  an on-screen keyboard will appear as shown below. The [Caps] button is used to switch between number/uppercase and lowercase number; the [Space] button is used to enter a space; the [Bkspace] button is used to delete the last character entered; press [OK] to confirm the entered data and exit the screen.

There may be a character limit depending on the content entered. Unavailable characters will be greyed out as you can see below.

In addition, you can select the numbers on the control panel for editing, press  to switch between numbers, lower and

upper case letters or characters and press  to delete the last character entered. Depending on the limits of the entered information, when changing the input method, selected characters are marked in green and unavailable characters are marked in grey.


**12.4 Data management**

In the main interface, select the button  to enter the data management screen, as you can see below. This screen shows all stored files. The user can use the search function (chapter 13.5) to select a file and edit it by deleting or modifying functions. In addition, the user can review the content of the stored files (chapter 12.6).

-  Go to the first page of the file list.
-  Go to the last page of the card list.
-  Go to the previous page of the card list.
-  Go to the next page of the card list.

**12.5 Searching for a File**


Select [Query] to open a data search window like the following. Enter the search conditions and press [Query] to obtain results. The [Clear] button is used to clear the inserted data. [Cond. and] and [Cond. or] indicate the search mode. The user can select either of them. If [Cond. and] is selected, the cards that match all the search data will be displayed. If [Cond. or] is selected, tokens matching one of the search data will be displayed.


Tip: It is recommended to use the [Cond. and] option, if there are many files saved. This way you can quickly find the file you are looking for.

### 12.6 Review of files

In the file management screen, hover over the desired file and select [Review] to open another window with the patient data as shown below. Users can modify the patient data. Select [Save] to save the changes, which will not be reversible.

To enter the review screen, select a file and press [Review]. A window similar to the sample interface will open as you can see below.

On the screen you can change the print mode by pressing . Press  to start printing.

**12.7 Date and time settings**

In the main interface, select the button  to open the date and time settings, as the following picture:

In this screen, you can change the options with the  buttons to edit the content.

**12.8 System Settings**

In the main interface, select the button  to open the system setting window, as in the following picture. In this dialogue window, select the [Default] button, to return to the default settings. Description of functions:


Option	Function	Description
Screen saver	None/30 sec/1 min/2 min/5 min/10 min	The screen saver will be activated when the selected time elapses without any action. "None" means that this function will not be used.
Backlight	30 sec/1 min/2 min/5 min/10 min/ Always ON	The backlight will turn off after the selected time has elapsed. "Always on means that it will not be switched off".
Light intensity	{90%}/[80%]/[70%]/[50%]	After selecting the light intensity, the display will light up with the set intensity.
Auto power off	1 min/3 min/5 min/10 min/15 min/30 min/60 min/none	The system will automatically switch off if no action is taken after the selected time has elapsed. "None" means that this function is not activated.
Low battery level	None/Once only/Always	The system will sound an alarm when the battery is about to run out.

Filter frequency	[50Hz/35Hz]/[50Hz/25Hz]/[60Hz/25Hz]/[60Hz/35Hz]	Setting the AC and EMG filter parameters.
Data input	Before]/[After]/[None] [Before]/[After]/[None] [English]/[Chinese]	Set when patient data is entered.
Language	[English]/[Chinese], etc.	Setting the system language.
K-B sound	On/Off	If you select "on", when you press the keys on the keyboard, the equipment will beep. If you select "off", it will not beep.
Demo Mode	Power On/Off	If you select "on", the system will operate in demo mode. Otherwise, it will operate in normal mode.

### 12.9 Sampling Configuration

Select the button  on the main interface to open the sampling configuration window. Select the [Default] button to return to the default configuration.


Description of the options:

Option	Function	Description
AC filter	On]/[Off] [On]/[Off] Function	AC filter on or off
EMG filter	[On]/[Off] [On]/[Off] Function	EMG filter on or off
DFT filter	[On]/[Off] [On]/[Off]	DFT filter on or off
Rhythm derivation	[I]/[III]/[III]/[avR]/[avL]/[avF]/[v1]/[v2]/[v3]/[v4]/[v5]/[v6]	Set the rhythm lead to use for printing in rhythm mode.
View	[3 leads]/[6 leads]/[6*2 leads]/[12 leads].	Set the ECG display view
Gain	[2.5mm/mV][5mm/mV][10mm/mV]/[20mm/mV][40mm/mV]	Set the gain
Rate	[5 mm/s]/[6.25mm/s]/[10mm/s] [12.5mm/s]/[25mm/s]/[50mm/s]	Set the rate. Automatic recording and rhythm recording are not supported with the 5mm/s, 6.25mm/s, 10mm/s, 12.5mm/s, 100mm/s options during printing.
Premature (%)	0~100	The system will use the value entered as the standard for diagnosing a premature beat.
Pause time (ms)	1200~3000	The system shall use the value entered as the standard for diagnosing a paused beat.
Tachycardia (bpm)	0~240	The system shall use the value entered as a standard for diagnosing tachycardia.
Bradycardia (bpm)	0~99	The system shall use the value entered as a standard for diagnosing tachycardia.

### 12.10 Print settings

Select the button on the main interface to open the print settings dialog box, as shown in the following figure:

Select the [Default] button to return to the default settings.  
The following table contains all objects, options and their explanation.

Option	Function	Description
Print mode	[Auto 6*2+1]/[Auto 6*2]/[Auto 4*3+1]/[Auto 3*4]/[Auto 2*6+1]/[Auto 2*6]/[Rhythm 7]/[Rhythm 6]/[Rhythm 5]/[Rhythm 4]/[Rhythm 3]/[Rhythm 2]/Manual	The selection shall be used as the default print mode.
Derivative gains	Intelligent/Current	The selection shall be used as the default bypass gains. "Smart" means that the system will adjust the gains automatically to adapt to the paper height during printing.
Auto strip	3 sec/4 sec/5 sec/6 sec/8 sec /10 sec/15 sec/20 sec/25 sec	The selected value will be the printing time duration for each line.
Rhythm strip	10 sec/15 sec/20 sec/25 sec/30 sec	In print mode, when rhythm 12, rhythm 10, rhythm 8 or rhythm 6 is selected, the system will use the selected value as the print duration for each row of waves.
Average QRS	[4*3+Mark]/[4*3]/[3*4+Mark]/[3*4]/[2*6+Mark]/[2*6]/[None]	In print mode, when "Rhythm" or "Auto" is selected, the system will use the selected format as the print duration for each row of waves
Self-diagnosis	All/Data/Conclusion/None	The self-test process contains 2 parts of data and conclusion. User can print one of the two.
Periodic	[for 1 min]/[for 2 min]/[for 3 min]/[for 5 min]/[for 10 min]/[for 20 min]/[for 30 min]/[for 60 min]/[off].	The system will periodically print ECGs at the selected time intervals. If the print mode option has been selected as "Manual", the system will print in "Auto 3x4+1". Otherwise, the system will print in the selected mode.

**Note:** Auto strip, Rhythm strip, QRS average, Auto-Diag and Periodicity options are only available if "Auto" or "Rhythm" print mode is selected.

### 12.11 Placement of electrodes

Select the button on the main interface to open the correct electrode placement image. Press any key to exit this interface.

## 12.12 Manufacturer information


Select the button  on the main interface to open the window with the equipment manufacturer information. This interface shows the name of the equipment, the version, the name of the manufacturer, copying rights and contact details of the manufacturer.

## 13. TROUBLESHOOTING

### 13.1 Automatic shutdown


- Check if the battery has run out.
- Check if the AC voltage is too high. The overvoltage control circuit causes the equipment to shut down automatically.
- Check if the AC interference is too high or if the plug screws are too tight. The automatic shutdown serves to protect the circuit.

### 13.2 AC Interference


- Is the earth wire of the device in good condition?
- Are the electrodes and cables properly connected?
- Are the electrodes and skin covered with sufficient gel?
- Is the metal bed installation adequate?
- Does the patient touch the wall or the metal bed?
- Is the patient touching other people?
- Are there any powerful electrical devices operating nearby? For example: X-ray or type B ultrasound equipment?

### 13.3 EMG Interference

- Check if the room is comfortable.
- Is the patient nervous?
- Is the bed too narrow?

### 13.4 Deviation from baseline

- A. Verify electrode adhesion and patient cable placement.
  - B. Check the connection between the patient cable and the electrodes.
  - C. Check the cleanliness of the electrodes and the patient's skin. Are the electrodes and gel covered with sufficient gel?
  - D. The patient should not move or hyperventilate.
  - E. Is the connection between the cable and the electrodes correct?
- We recommend using a filter if there is still interference after performing these checks.

### 13.5 Troubleshooting table

Problem	Cause	Recommended solution
Too much interference, wave moves.	<ol style="list-style-type: none"> <li>1. Incorrect earthing.</li> <li>2. Unstable cable connection.</li> <li>3. AC interference.</li> <li>4. Patient is nervous.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check the grounding wires.</li> <li>2. Place the patient in the correct position.</li> </ol>
Baseline is irregular.	<ol style="list-style-type: none"> <li>1. AC interference is too strong.</li> <li>2. Patient is nervous and there are strong EMG interferences.</li> </ol>	<ol style="list-style-type: none"> <li>1. Create a comfortable environment for the patient.</li> <li>2. If the bed is metal, replace the bed.</li> <li>3. The power cord and patient cable are not parallel or are too close together.</li> </ol>
Waveform is not regular, too large or too straight.	<ol style="list-style-type: none"> <li>1. Electrode conductivity is not good.</li> <li>2. The battery is flat.</li> <li>3. The contact between the electrodes and the skin is not adequate.</li> <li>4. The plug between the cable and the main unit is not properly connected.</li> <li>5. The contact between the patient cable and the electrodes is not correct.</li> </ol>	<ol style="list-style-type: none"> <li>1. Use high quality alcohol.</li> <li>2. Clean the electrodes and the area of the patient's skin that will come into contact with the electrodes.</li> <li>3. Charge the battery.</li> </ol>
Deviation from baseline.	<ol style="list-style-type: none"> <li>1. The battery has run out.</li> <li>2. Patient has moved.</li> </ol>	<ol style="list-style-type: none"> <li>1. Charge the battery.</li> <li>2. Instruct the patient to remain still.</li> </ol>
Waveform is not clear.	<ol style="list-style-type: none"> <li>1. The printhead is dirty.</li> <li>2. The paper is not suitable.</li> </ol>	<ol style="list-style-type: none"> <li>1. Clean the printer head with alcohol when the equipment is switched off, wait for the alcohol to volatilise before using the printer again.</li> <li>2. Use the appropriate thermal paper.</li> </ol>

## 14. MAINTENANCE

The customer is not allowed to open the device, to avoid electric shock. Any repairs or upgrades must be carried out by trained and authorised professionals. Maintenance must be carried out with original accessories. Please disconnect the power cable when switching off the device. If you are not going to use the device for a long time, please store it in a cool and dry place protected from sunlight and charge it once every three months.

**APPENDIX EMC INFORMATION**

**Guidance and manufacturer's declaration - electromagnetic emissions - for all EQUIPMENT AND SYSTEMS**

Guidance and manufacturer's declaration. Electromagnetic Emissions		
The electrocardiograph model ECG1200G (ECG) is intended for use in the electromagnetic environment specified below. The customer or the user of the electrocardiograph model ECG1200G (ECG) must use it in such an environment.		
Emissions	Compliance	Electromagnetic Environment: Guidance
RF Emissions CISPR11	Group 1	The electrocardiograph model ECG1200G (ECG) uses RF energy only for its internal operation. To this end, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions CISPR11	Class B	
IEC 61000-3-2 Harmonic Emissions	Class A	
Voltage fluctuations and flicker emissions IEC 61000-3-3	Complies	

**Guidance and manufacturer's declaration - electromagnetic immunity - for all EQUIPMENT AND SYSTEMS**

Guidance and manufacturer's declaration. Electromagnetic immunity			
The electrocardiograph model ECG1200G (ECG) is intended for use in the electromagnetic environment specified below. The customer or the user of the electrocardiograph model ECG1200G (ECG) should ensure that it is used in such an environment.			
Immunity Test	Test level 60601 EN/IEC	Compliance level	Electromagnetic environment: guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tiles. If the floor is covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast burst/transistor IEC 61000-4-4	±2 kV for power supply cables ±1 kV for input/output cables	±2 kV for power supply cables ±1 kV for input/output cables	Mains power quality should be typical of commercial or hospital environments.
Over current IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be typical of commercial or hospital environments.
Voltage dips, short interruptions and voltage variations on incoming power lines IEC 61000-4-11	<5% UT (>95% drop in UT) for 0.5 cycles 40% Ut (<60% drop in UT) for 5 cycles 70% Ut (<30% drop in UT) for 25 cycles <5% Ut (>95% drop in UT) for 5 sec.	<5% UT (>95% drop in UT) for 0.5 cycles Conforms 40% Ut (<60% drop in UT) for 5 cycles Conforms 70% Ut (<30% drop in UT) for 25 cycles Conforms <5% Ut (>95% drop in UT) for 5 sec.	Mains power quality should be typical of commercial or hospital environments. If the user of the equipment or system requires continuous operation during power interruption, it is recommended that the equipment or system be connected to a battery or uninterruptible power supply.
Power supply frequency (50/60 Hz) Magnetic field IEC 61000-4-8	3 A/m	3 A/m	Electrical frequency magnetic fields should be at levels typical of a commercial or hospital environment.
*Note: Ut represents the mains voltage prior to the application of the test level.			

Guidance and manufacturer's declaration. Electromagnetic Immunity
The electrocardiograph model ECG1200G (ECG) is intended for use in the electromagnetic environment specified below. The customer or the user of the electrocardiograph model ECG1200G (ECG) should assure that it is used in such an environment.

Immunity Test	Test level 60601 EN/IEC	Level of compliance	Electromagnetic Environment: Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz at 80 MHz	3	Portable and mobile RF communication equipment should not be used closer to any part of the ECG1200G (ECG) electrocardiograph, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5/\sqrt{V1}] \sqrt{P}$ $d = [3.5/E1] \sqrt{P}$ 80 MHz to 800 MHz $d = [7/E1] \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). The field strength of fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz at 2.5 GHz	3 V/m	

Note 1: At 80 MHz and 800 MHz, the upper frequency range applies.

Note 2: These guidelines may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a. Field strengths from fixed transmitters, such as radio base stations, telephones (cellular/wireless) and land mobile radios, amateur radio, AM and FM radio broadcasts cannot be predicted theoretically with accuracy. To assess the electromagnetic environment caused by fixed RF transmitters, an electromagnetic analysis of the site should be considered. If the measured field strength at the location of the equipment exceeds the applicable RF compatibility level mentioned above, it is recommended to observe the behaviour of the monitor to verify that it is operating normally. If abnormal behaviour is observed, additional measures, such as re-orientation or relocation of the device, may be necessary.
- b. In the frequency range 150 KHz to 80 MHz, field strengths should be less than 3V/m.

**Recommended separation distances between portable and mobile RF communications equipment and the equipment or system - for non-life-supporting equipment or systems.**

Recommended separation distances between portable and mobile RF communications equipment and Model ECG1200G (ECG)			
The Model ECG1200G (ECG) electrocardiograph is intended for use in an electromagnetic environment in which radio frequency interference is controlled. The customer or the user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile radio frequency communications equipment (transmitters) and the unit recommended below, depending on the maximum output power of the communications equipment.			
Maximum transmitter power [W] 0.01	Distancia de separación d [m] según la frecuencia del transmisor		
	150 kHz a 80 MHz $d = [3.5/\sqrt{V1}] \sqrt{P}$	80 MHz a 800 MHz $d = [3.5/\sqrt{V1}] \sqrt{P}$	800 MHz a 2,5 GHz $d = [7/E1] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	0.73
100	12	12	23
For transmitters whose maximum output power is not included in this table, the recommended separation distance d in metres (m) can be calculated from the equation applicable to the frequency of the transmitter, where P represents the maximum output power in watts (W) according to the transmitter manufacturer.			
Note 1: at 80 MHz and 800 MHz the upper frequency range applies.			
Note 2: These guidelines may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

 CONTEC MEDICAL SYSTEMS CO.,LTD  
 No.112 Qinhuang West Street, Economic & Technical Development Zone, Qinhuangdao, Hebei Province, PEOPLE'S REPUBLIC OF CHINA

 Legal representative inside EU/Representante legal en la UE: Shanghai International Holding Corp. GmbH (Europe)  
 Eiffestrasse 80, 20537 Hamburg, Germany

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Imported by/Importado por: Grupo R. Queraltó S.A.