ELECTROCARDIÓGRAFO


ECG300G

Manual de instrucciones

ES

Instruction manual

EN

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
ECG300G

NOTICE

Notice to the user and/or patient: any serious incident which has occurred in connection with the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Content

Chapter 1. Main Technical Specifications.....	1
Chapter 2. Safety Notice.....	2
Chapter 3. Maintenance Standards.....	2
Chapter 4. Instrument Characteristics.....	3
Chapter 5. ECG300G Panel Sketch Map.....	4
Chapter 6. Care Before Operation.....	6
Chapter 7. Preparatory Work Before Instrumentation.....	6
Chapter 8. Caution During Operation.....	6
Chapter 9. Electrode Placement.....	7
Chapter 10. Grounding and Connecting the Power Supply to the Instrument.....	8
Chapter 11. Caution for Battery Operation	9
Chapter 12. Control Panel and Key Instruction.....	9
Chapter 13. IR Troubleshooting.....	10

Chapter 1. Main Technical Specifications

1.1 Normal working environment

Operation

- a) Ambient temperature: +5°C~+35°C
- b) Relative Humidity: ≤80%.
- c) Power supply: AC:220V, 50Hz(110V, 60 Hz)
DC: 7.4V, 3700 mAh rechargeable lithium battery
- d) Atmospheric Pressure: 86kPa~106kPa

Storage and Transport

- a) Ambient temperature: -10°C~55°C
- b) Relative Humidity: ≤95%.
- c) Atmospheric pressure: 50kPa~106kPa

1.2 Input: Float and defibrillation protection

- 1.3 terminals: Standard 12 terminals
- 1.4 Patient leakage current: <10μA
- 1.5 Input Obstruction: ≥50MΩ
- 1.6 Frequency response: 0.05Hz~150Hz (-3dB)
- 1.7 Time constant: Time constant >3.2s
- 1.8 CMRR: >60dB, >100dB (Add filter)
- 1.9 EMG interference filter: 35Hz(-3dB)
- 1.10 Recording method: thermal printing system
- 1.11 Specification of recording paper: 80mm(W)×20m(L) High-speed thermal paper

1.12 Paper speed:

- Automatic recording: 25mm/s, 50mm/s, error:±5%.
- Pace registration: 25mm/s, 50mm/s, error:±5%.
- Manual registration: 5mm/s, 10mm/s, 12.5mm/s, 25mm/s, 50mm/s, error:±5%.

1.13 Sensitivity selection: 5, 10, 20mm/mV, error:±5%. Standard sensitivity is 10mm/mV±0.2mm/mV.

1.14 Auto recording: Record following recording format and auto mode, auto lead change, auto measurement and analysis.

- 1.15 Pace recording: The process records following pace format, pace mode, auto measurement and analysis.
- 1.16 Manual recording: Record following recording format, manual lead change.
- 1.17 Measurement Parameters: HR, P-R interval, P Duration, QRS Duration, QRS Duration, T Duration, Q-T interval, Q-Tc, P-Axis, QRS Axis, T-Axis, R(V5), S(V1), R(V5)+S(V1)
- 1.18 Product Safety Type: Class I, Type CF, defibrillation circuit and pacing protection circuit are present.
- 1.19 Bias Voltage Supported: $\pm 300\text{mV}$
- 1.20 Noise level: $\leq 15\mu\text{Vp-p}$
- 1.21 Fuse specification: 2pcs $\phi 5 \times 20\text{mm}$ AC range; T250mA/250V(Power supply:220V)
- 1.22 Size: 315mm(length) \times 215mm(width) \times 77mm(height)
- 1.23 Net weight: 2.25Kg

Chapter 2. Safety Notice

- 2.1 The electrical power supply must be properly grounded before operation.
- 2.2 If the grounded cable is not integrated, the device has to be operated with a built-in power supply.
- 2.3 Please pull out the power supply plug before changing the fuse.
- 2.4 This device must be operated and stored by professional personnel.
- 2.5 The operator must read this user's manual carefully before operating, and operate the device strictly in accordance with the operation rules.
- 2.6 The design of this device has been made based on a mature consideration of safety, but the operator should never neglect the condition of the device and the patient's situation.
- 2.7 Please switch off the equipment and remove the plug from the power supply before cleaning and disinfection.
- 2.8 Please do not operate this device in environments containing flammable anaesthesia gas.
- 2.9 If the device is used together with a cardiac defibrillator or other electrical stimulation devices at the same time, please select Ag/AgCl chloride chest electrode and ECG lead with the preventive defibrillation function. To prevent the metal electrode from burning the patient's skin, the disposable chest electrode should be used if the defibrillation duration exceeds 55 seconds. It is best not to use this device in conjunction with another electrical stimulation device at the same time. If there is no other choice but to use them at the same time, professional technical guidance is professional technical guidance is essential at the scene.
- 2.10 When other devices are connected to this ECG instrument, they must be Type I devices in accordance with IEC60601-1. As the total amount of leakage current may harm patients, the monitoring leakage current is carried by the connected devices.
- 2.11 The following descriptions relate to special attentions in ECG measurement and interpretation.
 - 1) The P wave and Q wave are not always reliable in the Muscle Intensive Artifact or Alternating Current Interference file. The same is true for the ST segment and T wave.
 - 2) Unclear terminals and windings of the S-wave and T-wave can lead to tolerance in the measurement process.
 - 3) In the archive, the R wave has been left out due to the low voltage of the QRS wave or any leads that have already fallen out. The measured heartbeat may deviate greatly from the correct value.
 - 4) Axis calculation and QRS border line identification are not always reliable in the low QRS wave voltage file.
 - 5) Occasionally, frequent ventricular premature complexes may be identified as beat dominant rhythm.
 - 6) Versatile arrhythmia fusion may result in unreliable measurement due to the difficulty in distinguishing the P wave in such a situation.
 - 7) ECG300G has been designed to perform ECG trace interpretation immediately after measurement. It is precisely this interpretation that does not report all possible heart problems and may sometimes not comply with the physician's diagnosis. Therefore, the final conclusion regarding each patient is up to the physician based on the patient's symptoms, the interpretation of the ECG300G unit and other examinations.**


Chapter 3. Maintenance Standards

- 3.1 In normal use condition according to the user manual and operation notices, if this instrument has a problem please contact our customer service department. Our company has the sales record and customer files for each instrument. The customer has one year warranty service from the day of shipment according to the following time and conditions. In order to be able to supply and service our instruments, we have a one year warranty service from the day of shipment according to the following time and conditions.
- 3.2 Our company can adopt ways of instruction, sending mail to the company by courier, visiting customers' companies, etc. to carry out the maintenance promise.
- 3.3 Even in the period of free maintenance, there is a charge for repairs in the following files:
 - 3.3.1 Failure or damage caused by misuse by not operating in accordance with the user's manual and operation notices.
 - 3.3.2 Failure or damage caused by accidentally dropping when moving after purchase.
 - 3.3.3 Failure or damage caused by preparation, reconstruction, decomposition, etc. outside our company.
 - 3.3.4 Failure or damage caused by natural disasters such as fire, flood, earthquake, etc.
 - 3.3.5 Failure or damage caused by unsuitable thermal recording paper.

- 3.4 The free maintenance period for spare parts is half a year. Power cable, recording paper, operation manual and packing material are excluded.
- 3.5 Our company is not liable for failures of other connecting instruments caused by failures of this device directly or indirectly.
- 3.6 The free maintenance service will be cancelled if the protection label has been destroyed.
- 3.7 For charge maintenance beyond the warranty period, our company advises to continue using "Maintenance Contract Standard". Please consult our customer service department for specific situations.




Chapter 4. Characteristics of the Instrument

- 4.1 Registration system: Matrix-thermal (8 dots/mm). You do not have to adjust anything. Frequency response is 150Hz maximum.
- 4.2 The device can record waveform and three-channel ECG observation, clear and accurate in real time, continuously. The observation includes: sign of power supply, sensitivity, paper speed, filtering status, etc..
- 4.3 In automatic mode, press the button once, start the recording process, which can improve your working efficiency.
- 4.4 On-screen keyboard control, more convenient for operation. The touch screen displays the working status, more clear for observation.
- 4.5 Security Class: Class I, Type CF.
- 4.6 Power supply includes both AC and DC power supply. The device includes a built-in rechargeable lithium battery.
- 4.7 This instrument is able to record 150 pieces of ECG waveform and print for 90 minutes continuously under the best DC status.
- 4.8 This instrument is able to store more than 1000 pieces of patient data, more convenient for data review and statistics.
- 4.9 The figure of the whole device is elegant and smooth.
- 4.10 According to the degree of defense of deleterious fluid, this device is common type.
- 4.11 According to the degree of safety used under the condition with anaesthetic gas mixed with air (or oxygen, nitrous oxide), this device belongs to device which cannot be used under the condition with flammable anaesthetic gas, mixed with air (or oxygen, nitrous oxide).
- 4.12 Digital signal processor for effective inhibition of baseline drift, interference, etc.
- 4.13 The instrument has the function with regular automatic measurement of ECG waveform parameters, self-analysis and self-diagnosis. It will help to reduce the abundance of doctors and improve the working effect.
- 4.14 According to the kind of working mode, this device is non-continuous working type.
- 4.12 Explanation of some symbols on this device:

Symbol	Description
~AC	AC working mode
OFF	Power supply disconnected
ON	Connected power supply
	Equipotential point
	Sites to be aware of, please refer to the user manual.
	The device type is CF, with defibrillation protection function.
	RS232 Connector
	USB Connector
 PATIENT	Cable connector

Chapter 5. ECG300G Panel Sketch Map

5.1 Outline map and name of components

5.2 Definition of the Key

-  — Function key: On/Off
-  — Function key: Positive adjustment
-  — Function key: Paper speed adjustment
-  — Function key: Filter function selection
-  — Function key: Menu
-  — Function key: Change print mode
-  — Function key: 1mV marker
-  — Function key: Print
-  — Function key: Enter
-  — Function key: Up
-  — Function key: Down
-  — Function key: Left
-  — Function key: Right

5.3 Definition of Indicator

The indicator turns green when AC power is supplied, and when the indicator turns green and red at the same time it is recharging.

Indicator for instrument when switched on.

Chapter 6. Care Before You Operate

- 6.1 You are required to read this operation manual thoroughly before operating to ensure safe and effective operation of this equipment.
- 6.2.1 No high voltage cables, radial X instruments, ultrasound instruments, electrotherapeutic instruments, etc., should be located around the ECG device.
- 6.2.2 Do not use or reserve the instrument in any place where the atmospheric pressure is too high, where the temperature and humidity exceed the common standard, where the ventilation is not satisfactory and there is too much dust or where there is gas containing salt and alkaline and chemical medicine.
- 6.3 The instrument should be placed on a flat surface. Pick it up and place it gently when moving it. Avoid excessive vibration and shock.
- 6.4 The value of AC current frequency and voltage must be in accordance with requirements and sufficient current capacity must be guaranteed.
- 6.5 Please place the device in a place where the grounding is easy to be done. Do not connect patients, do not connect patient connection cables with other conductors including earth, or beds which can be well conducted to earth.

Chapter. 7 Preparatory Work before Instrumentation

- 7.1 Check that the instrument is properly grounded and that the cable connections are secure.
- 7.2 Ensure that all electrodes connected directly to patients are in proper and firm condition.
- 7.3 Check the output voltage when selecting DC UPS.
- 7.4 Pack the gel separately, avoiding short-circuit caused by mutual touch between chest electrodes.
- 7.5 The AC power cable cannot be linked with the ECG cable.

Chapter 8. Caution during Operation

- 8.1 Pay attention to the condition of the patient and the instrument at all times.
- 8.2 Both the patient and the instrument may be connected only with ECG leads.
- 8.3 Maintain close observation of the patient and the instrument, ensuring that the patient and the instrument do not
- 8.3 Keep close observation of the patient and the instrument, ensuring that they do not move during the operation.
- 8.4 Turn off the instrument after use.
- 8.5 Disconnect the power and remove the ECG leads gently, without applying excessive force.
- 8.6 Correctly and properly store the instrument and spare parts for operation the next time.
- 8.7 Paper Loading


8.7.1 The dimension of the high-speed thermal printing paper in this instrument is: 80mm(W)×20m(L).

8.7.2 Open the paper cabinet cover, pull out the paper shaft, and install the registration paper tightly until the figure enters the correct internal position.

8.7.3 Close the paper cabinet cover. It is recommended to leave 2cm of registration paper outside.

Chapter 9. Registration Paper Instruction

9.1 The message “No Paper” will be displayed on the LCD when any recording paper is in shortage.

9.2 Specific high sensitivity paper is recommended for a clearly permanent print and may damage the printing mechanism. Please consult the distributor or the manufacturer for more details on how to purchase the paper.

9.3 Failure of the recording paper may be affected by high temperature, bad humidity or direct sunlight. For long term storage, the recording paper should be placed in a dry, dark and cool area.

9.4 Substance that may cause staining of the recording paper: Gel, glue and wet composite paper including its organic solvent.


9.5 Substance that can cause waves to fade:

File folders made of soft PVC material, plastic, etc; eraser and magnetic tape containing plasticiser. Fluorescence and stamping ink.

Chapter 10. Electrode Placement

Note: Set up the chest electrode first and then the limb electrode.

10.1 Chest Electrode (See Figure 4)

At

ow:

V1: Inner fourth costal to the right edge of the sternum.

V2: Inner fourth costal to the left edge of the sternum.

V3: Midway between V2 and V4.

V4: Fifth internal costal at the left mid clavicular line.

V5: Left anterior axillary line at the horizontal level of V4.


V6: Left middle axillary line at the horizontal level of V4.

Thoroughly clean the skin where the chest electrodes will be attached with alcohol and then apply ECG cream to the site, over 25mm in diameter and also to the edge of the chest electrodes. Press and attach the electrodes to the positions from V1-V6.

Keep in mind that it is not allowed to let the electrodes come in contact with each other or for the cream to overlap from one position to another.

10.2 Limb Electrode

Electrodes should be placed on the soft skin of the hands and feet. Clean all limb electrodes and the positions around where the limb electrodes are to be attached with alcohol before applying ECG cream, then firmly snap the electrodes into the positions.

Caution: Secure the handle of the ECG cable plug tightly after inserting it into the instrument.

10.3 Checklist for Electrode and ECG Cable Connection

Location of Electrode	Electrode Code	No. of Slot
Right Arm	RA/R	9
Left Arm	LA/L	10
Left Leg	LL/F	11
Right Leg	RL/N	14
Chest 1	V1/C1	12
Chest 2	V2/C2	1
Chest 3	V3/C3	2
Chest 4	V4/C4	3
Chest 5	V5/C5	4
Chest 6	V6/C6	5

Chapter 11. Grounding and Connecting the Power Supply to the Instrument Instrument

Make sure that the instrument status is off and then put the ground phase on the instrument through a 3-prong outlet. An earth wire can be used to connect the ground terminal of the instrument. Do not use other pipes. The proper ground phase can guarantee safety and prevent AC power and electromagnetic wave interference.

Chapter 12. Battery Operation Caution

12.1 This instrument is designed with a built-in, maintenance-free, sealed, rechargeable lithium battery. It has an automatic charge and discharge monitoring system. The battery shall be automatically recharged when the instrument is connected to the AC mains power supply.

12.2 The device is capable of printing continuously for 90 minutes and working 4 hours without printing after fully charging the battery. When it is working, the LCD displays the battery status signal in 5 degrees. When the battery lacks electric power, the instrument will automatically shut down to avoid damage to the battery.

12.3 The battery should be recharged in time after a period of strenuous use. For long term storage, it is required to recharge the battery once in every 3 months. The life of the battery can be extended if this is done.

12.4 Seven statuses of the battery power supply are displayed on the LCD as follows:

No.	Mark	Description
b	■ ■ ■	Unknown status, normally displayed when just logged in to the instrument within 1 minute.
c		Using an alternating current power source
d		Using battery, and full power
e		Using battery, volume: 3/4
f		Using battery, volume: 1/2
g		Using battery, volume: 1/4
<p>Note: When charging the battery, the icon changes from f to c. Using battery, but low power, it is suggested to recharge the battery or use the AC power supply.</p>		

12.5 When the battery cannot be recharged or cannot work for more than 10 minutes after being fully charged, please replace the battery with another one.


WARNING!!!

- Do not connect either the polar “+” or the polar “-” of the battery directly to the wire, otherwise a fire hazard will occur.
- Danger of explosion is possible if the instrument is stored near a burning lighted area.
- Do not open or disassemble the battery.

Chapter 13. Control Panel and Key Instruction

13.1 Main Interface


It is displayed as follows:

Electrical Power Status: Please refer to 12.4
Keyboard:


	Enter the sampling interface. When the instrument is in the on state, it will start this operation automatically
	Enter the File management interface. Query, modify or delete file information.
	See the sketch map for electrode placement information.
	Date and Time settings.
	System Configurations.
	Sampling configurations.
	Analysis parameter settings, settings for each parameter used for automatic analysis.
	Print settings, print mode settings, style and content.
	About us, view information about our company and software version.

13.2 Sampling Interface

Select the button  in the main menu or choose the shortcut key  to enter the sampling interface.

Attention: Due to the “configuration”, Patient information can be entered before the sampling signal, appearing with the option: entering file information.

Multi-form terminal waveforms are displayed, including 3 terminals per screen, 6 terminals per screen and 12 terminals per screen. The 12-terminal style sampling interface can be displayed as follows: 1.

1. Print Style
2. Terminal Off
3. Status
4. Magnification
5. Speed
6. Filter
7. Current heartbeat

Stop Sampling: Press the key  on the keypad to stop the sampling process and return to the main interface.

Change of Terminal: You can press the  /  key to display the other terminals. The terminals are used for manual printing.

Change Terminal Style: Press the  /  key to change the terminal style between 3 terminals per display, 6 terminals per display and 12 terminals per display.

Terminal Off: In Demo mode, “DEMO ECG” is displayed. In sampling mode, the Terminal Off information is displayed.

Change Registration Style: Press the key  to change the recording mode between auto 3x4+1, auto 3x4, auto 2x6+1, auto 2x6, auto 2x6, auto 4x3, rhythm 4, rhythm 3, rhythm 2 and manual mode.

Adjust Magnification (Sensitivity Selection): The sensitivity is selected by pressing the key . User can select 5mm/mV, 10mm/mV or 20mm/mV.

Adjust the speed: Use the speed adjustment key  to change the speed: 5mm/s, 10mm/s, 12.5mm/s, 25mm/s, 50mm/s. Auto logging and rhythm logging cannot support 5mm/s, 10mm/s, 12.5mm/s.

Change the filter: Use the filter selection key  to switch between no-filter, AC, EMG, DFT, AC+EMG, AC+DFT, EMG+DFT, AC+EMG+DFT.

Display calibration: Use the “1mV”  key to display the 1mV marker on the screen.

Print/Finish printing: Use the “print”  key on the keyboard and then a print operation is started or stopped.


Automatic Mode: After starting printing, the system will print and store 12 synchronous waveforms automatically. The length is decided by the setting in the print option. And also according to the setting item, the automatically analysed data and results will be printed and the system can end the printing automatically.

Manual Mode: After starting the printout, the user must print the actual waveform by switching between terminal groups. It means that the ECG waveform printed under manual mode is not synchronous and cannot be stored. And the user must stop the printout by pressing the key again.

During the printing process, the printing status includes:

Display content	Observation
Processing...	in the printing process
Waiting...	in the process of finishing printing
No paper	Lack of paper. The user has to restart the operation after filling in papers
Out of print lead time	The connection between the system and the printing sub-system is broken.
Out of ECG waiting time	The connection between the system and the sampling sub-system has been broken.
Low electrical power	With low power, the system cannot start the print job.

On the interface, press the button **SET** on the keypad and the system will enter the shortcut configuration interface.

Select [OK] button and the system will apply the new settings and return to the sampling interface. Select [Cancel] and the system will return to the sampling interface without applying the new settings.

Each of the option functions has been shown in the following table:

Item	Optional Content	Observation
AC Filter	(On)/(Off)	Setting whether to use the AC Filter or not
EMG Filter	(On)/(Off)	Setting whether to use the EMG Filter or not
DFT Filter	(On)/(Off)	Setting whether to use the DFT Filter or not
Rhythm Terminal	Any of the 12 terminals	Setting up the rhythm terminal to print ECGs in rhythm print mode
Show Style	(3 terminals)/(6 terminals)/(12 terminals)	Wave configuration, show style
Show Augmentation	(5mm/mV)/(10mm/mV)/(20mm/mV)	Wave configuration, show magnification
Show Speed	(5mm/s)/(10mm/s)/(12.5mm/s)/(25mm/s)/(50mm/s)	Waveform configuration shows speed. Automatic recording and rhythm recording cannot support 5mm/s, 10mm/s and 12.5mm/s.

13.3 Entering File Information

According to the different configuration items (see 13.8), the user can enter the patient file before or after sampling and can also enter a blank file. The input window is as follows:


Select any input window. Press the key **SET**, the "On-Screen Keyboard" will appear as follows. The function of the [Caps] button on the "On Screen Keyboard" is to change the key from lower number and file to upper score and file. Press [OK] to confirm the input and exit this interface.

There may be a character limit depending on the content entered. And then the limited character(s) will be greyed out and unavailable.

13.4 Archive Management

In the main interface, select the button  and you can enter the archive management interface, as shown:

Sum of current files / Sum of total files

Currently selected file

List of archives

Display menu advanced

Page roll buttons

Delete the selected file selected

Exit this menu

Display the revision dialog box

Date and time	ID	Name	Sex
2008-06-11 08:30:53	0001	GaoXiaohua	F
2008-06-10 15:31:15	0002	SongGuang	M
2008-06-10 15:12:39	0003	LiCui	F
2008-06-10 14:51:26	0004		
2008-06-10 14:27:21	0005		
2008-06-10 14:25:03	0006	ZhangXiaohua	M
2008-06-04 14:08:49	0007		
2008-06-04 11:57:12	0008	QinXiaohua	F
2008-06-04 09:33:02	0009		
2008-06-04 09:14:47	0010	CaiLinlin	M

This interface displays all stored files. Users can use the search function (refer to 13.5 file inquisition) to select the required file, and edit any file by delete or modify operation. In addition, the user can review the stored file information. (refer to 13.6 file review).

- Go to the first page of the archive list.
- Go to the last page of the file list.
- Go to the previous page of the file list.
- Go to the next page of the file list

13.5 Archive Query

Select [Adv-Opr] in the archive management and then, the following window appears


Advanced menu

Date and time	ID	Name	Sex
2008-06-13 10:02:58			
2008-06-13 09:55:01			
2008-06-11 15:24:05			F
2008-06-11 08:30:53	0001	GaoXiaohua	F
List All	1:15	0002	SongGuang
Query	2:39	0003	LiCui
Delete All	7:21	0005	
Return	5:03	0006	ZhangXiaohua
	3:49	0007	

Select [Query] to launch a File Query dialog box as shown below.

Enter search condition and select the [Select] button and the user can have the desired result. The function of [Clear] is to clear the entered query condition.


[Cond. And] and [Cond. O] indicate the match mode of the search condition, and the user can select either one. If [Cond. And] has been selected, the search result will match all the conditions entered simultaneously. If [Cond. Or] has been selected, the search result will display the factors that match any of the conditions.


Tip: If you are looking for a large number of patient files, enter the confirmed search conditions and select [Cond. And] and you will be able to find the desired patient files immediately.

13.6 File Review

In the File management interface, after moving the focus to the correct patient file being reviewed, select [Review] to launch the following dialog box where the patient file information is displayed. Users can modify here, select [Save]. The change, which is irreversible, will be saved.

Be sure to make the correct selection. Select [Review] to start the next review interface, which is similar to the sampling interface.

In this interface, the user can adjust the time segment of the waveform by means of the ◀ and ▶ button and also the speed and gain. (refer to 13.2 sampling interface). In this state, user can print by pressing . By pressing **SET**, the system shall enter the shortcut configuration interfaces as follows:

Select (OK) button, the system will apply the new settings and return to the review interface. Select (Cancel) button, the system will return to the review interface without applying the new settings.


Each of the option functions is shown as follows:

Item	Optional Content	Optional Content
Print mod	(Auto 3x4)/(Auto 2x6)/(Auto 4x3) and any mode that fits the current file.Rhythm terminal	Print mode settings
Rhythm terminal	Any 12 terminals	Configure the rhythm terminal to print ECGs in rhythm print mode.
Show style	(3 terminals)/(6 terminals)/(12 terminals)	Wave configuration, show style
Show willingness	(5mm/mV)/(10MM/Mv)/(20mm/mV)	Wave configuration, show magnification
Show speed	(5mm/s)/(10mm/s)/(12.5mm/s)/(25m/s)/(50mm/s)	Wave configuration, display speed. Automatic recording and rate recording cannot support 5mm/s, 10mm/s, 12.5 mm/s.

13.7 Date and Time Settings

In the main interface, select button  to open the following dialogue window where the Date and Time settings are displayed.


In this interface, users can select the key  and  to switch between all items, using  and  to adjust the content of options.

13.8 System Settings

In the main interface, select the button  to open the following dialogue window with system settings.

In this dialog box, select the (Default) button, the system settings will return to default parameters. Each of the options functions as follows:

Item	Item	Observation
Screen Protector	None / 30 seconds / 1 minute / 2 minutes / 5 minutes / 10 minutes	The screensaver will be archived after the selected period of time. "Nothing" means that this function will not be used.
Backlight	30 seconds / 1 minute / 2 minutes / 5 minutes / 10 minutes / Always On	The backlight will be switched off after the selected period of time. "Always On" means that the backlight will never be switched off.
Auto Power Off	None / Only once / Always	The system will be shut down if there is no operation after a selected period of time. Nothing means that the function is not effective.
Low Power	Before / After / Nothing	Which alarm regime the system will adopt when the battery power is about to run out.
Input Information	[50Hz/35Hz]/[50Hz/25Hz]/[60Hz/25Hz]/[60Hz/35Hz]Filter	Set the time and date when you enter the file information.
Filter Frequency	[English]/[Chinese], etc.Language	Set the AC Filter and EMG Filter parameter
Sound K-B	On / Off	When you press the key on the keyboard, the instrument will produce a sound if "On" has been selected. Otherwise, no sound will be produced.
Demo Mode	On / Off	The system will operate under the demo version, if "On" has been selected. Otherwise, it will operate under the normal version.

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Select the [Default] button and the sampling settings will return to the default settings. Each of the option functions is shown as follows:

Item	Optional Content	Observation
AC Filter	On/Off	Default setting using AC Filter or not
EMG Filter	On/Off	Default setting using EMG Filter or not
DFT Filter	On/Off	Default setting using DFT Filter or not
Show Style	[3 terminals]/[6terminals]/[[12 terminals]	Default settings show style
Show magnification	[5mm/mV]/[10mm/mV]/[20mm/mV]	Default setting show magnification
Show feed	[5mm/s]/[10mm/s]/[12.5mm/s]/[25mm/s]/[50mm/s]	Default settings show speed. Auto log and pace log cannot support 5mm/s, 10mm/s and 12.5mm/s.

13.10 Analyse Parameter Configurations

Select the button  in the main interface to launch the following dialogue window in which the parameter settings are analysed.

The settings will affect the real time analysis diagnostic track, file review and report printing during the sampling process.

Rhythm Lead:

Premature(%):

Pause Time(ms):

Tachycardia(bpm):


Bradycardia(bpm):

Select the [Default] button, the system settings will revert to the default settings.
Refer to the following:

Item	Observation
Pace Terminal	Set the rhythm terminal to analyse heart rate and print ECGs in rhythm print mode.
Premature	The system will use the value entered as a standard for considering the premature beat.
Pause time	The system will use the value entered as a standard for considering the pause beat.
Tachycardia	The system will use the value entered as a standard for considering tachycardia.
Bradycardia	The system will use the value entered as a standard for considering bradycardia.

13.11 Print Settings

Select the button  on the main interface to launch the print settings dialog box, as shown:


Select the [Default] button, the print settings will return to the default settings.
In this interface, the auto mode option can be effective only when you select “auto” in [Print Mode].

Item	Optional Content	Observation
Print mode	[Auto 3x4+1]/[Auto 3x4]/[Auto 2x6+1]/[Auto 2x6]/[Auto 4x3]/[Rythm2]/[Rythm3]/[Rythm4]/[Manual]	The selection will be used as the default print mode.
Cable magnification	Smart/Current	The selection will be used as the default Cable Augmentation. "Smart" means that the system will adjust the wire rise automatically to match the paper height during printing. "Actual" means that the system will use the wire magnification on the screen during printing.
Self act of undressing	3seconds/4 seconds/5 seconds/6 seconds/8 seconds/10 seconds/15 seconds/20 seconds/25 seconds	The selection will be used as the duration of printing time per default.
Stripping rhythm	10 seconds / 15 seconds / 20 seconds / 25 seconds / 30 seconds	With the print mode set to "Rhythm 2", "Rhythm 3" or "Rhythm 4", the system will use the selected option as the duration of the print time of the stripping rhythm.
QRS average	[2x6]/[2x6+Marca]/[3x4]/[3x4+Marca]/[4x3]/[4x3+Marca]/[Nada]	With the print mode is "Rhythm" or "Auto" or "Rhythm 4", the system will use the selected option as the print time duration of the stripping rhythm.
Self-diagnosis	All / Fact / Conclusion / Nothing Periodic	The self-diagnostic process contains 2 parts of data and conclusion. User can print one of the two as desired.
Periodic	[for 1 min]/[for 2 min]/[for 3 min]/[for 5 min]/[for 10 min]/[for 20 min]/[for 30 min]/[for 60 min]/[Off].	The system will periodically print ECGs at the selected time interval. If the print mode option has been selected as "Manual", the system will print in "Auto 3x4+1". Otherwise, the system will print in the selected mode.

Note: Auto Act Stripping, Rhythm Stripping, QRS Averaging, Autodiagnostic, Periodic are available when the "Auto" or "Rhythm" print mode has been selected.

13.12 Check Electrode Placement

Select the button  on the main interface to start the following power line placement illustration interface.

Select any key to exit this interface

13.13 About us

Select the button on the main interface and you will be able to start the next interface of information related to this instrument. This interface shows the instrument name, version, company name, copy right and company contact details.

Contec Medical Systems Co., Ltd.
No. 112 Qinhuang West Street, Economic &
Technical Development Zone, Qinhuangdao,
Hebei Province, PRC

Legal representative inside EU/
Representante legal en la UE

Shanghai International Holding Corp.
GmbH (Europe) Eiffestrasse 80,20537
Hamburg, Germany

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Imported by/Importado por: Grupo R. Queraltó S.A.