# ESPIRÓMETRO

## SP10

Manual de instrucciones

ES

Instruction manual

EN

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SP10

**NOTICE**

**Notice to the user and/or patient: any serious incident which has occurred in connection with the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.**

## Users Manual

### Warning:

- For accuracy it is recommended that the SPIROMETER should not be tested on the same examinee more than 5 times.
- The examinee should blow out all the air during the test, do not change air or cough.
- Do not use the device in a lower temperature environment.
- Automatic shut-off when not in use for one minute.
- Please refer to the correlative literature on clinical restrictions and caution.
- This device is not intended for treatment.

The User Manual is published by our company. All rights reserved.

## Content

Chapter 1. Safety.....	2
1.1 Instructions for safe operation.....	2
1.2 Warning .....	2
1.3 Attention.....	2
Chapter 2. Overview.....	2
2.1 Characteristics .....	3
2.2 Major applications and outreach.....	3
2.3 Environmental requirements.....	3
Chapter 3. Principle.....	3
Chapter 4. Technical specifications.....	3
4.1 Main performance.....	3
4.2 Main parameters.....	3
Chapter 5. Installation.....	4
5.1 Front panel view.....	4
5.2 Installation.....	4
5.3 Accesories.....	4
Chapter 6. Operation Guide.....	4
6.1 Method of application.....	4
6.1.1 Installation.....	4
6.1.2 Measurement.....	4
6.1.3 Parameter interface.....	5
6.1.4 Menu operations.....	5
6.1.5 Repeated measurement.....	10
6.1.6 Upload.....	10
6.1.7 Load data.....	10
6.2 Attention.....	11
Chapter 7. Maintenance, transport and storage.....	11
7.1 Cleaning and disinfection.....	11
7.2 Maintenance.....	11
7.3 Transport and storage.....	11
Chapter 8. Troubleshooting.....	11
Chapter 9. Key symbol.....	12
Chapter 10. Entering parameters.....	13

## Chapter 1. Safety

### 1.1 Instructions for Safe Operations

- Check the main unit and all accessories periodically to make sure that there is no visible damage that may affect patient's safety and monitoring performance. It is recommended that the device should be inspected weekly at least. When there is obvious damage, stop using the device.
- All maintenance must be performed by qualified service engineers ONLY. Users are not permitted to maintain it by themselves.
- The SPIROMETER cannot be used together with devices not specified in User Manual. Only the accessory that is appointed or recommendatory by manufacture can be used with this device.
- This product has been calibrated before leaving factory.

### 1.2 Warning

- Explosive hazard- DO NOT use the SPIROMETER in the environment with tinder such as anesthetic.
- Please check the packing before use to make sure the device and accessories are totally in accordance with the packing list, or else the device may have the possibility of working abnormally.
- Don't use the device environment with strong electromagnetic interference, direct breeze source, cold source and hot source.
- Portable or mobile RF equipment with strong electromagnetic interference may influence the accuracy of this device.
- Improper disposal of device and its accessories and packing (include mouthpiece, plastic bags, foams and paper boxes) may cause environment pollution, please follow the local laws and regulations.
- Please choose the accessories which are appointed or recommended by the manufacturer for avoiding device damage.
- Don't use the device with the turbine or the same kind product.
- DO NOT use the device when it is under charging state.
- The red and green indicators are all highlight in charging state, the red indicator goes out when the charge has finished.

### 1.3 Attention

- Keep the SPIROMETER away from dust, vibration, corrosive substances, tinder, high temperature and moisture.
- If the SPIROMETER gets wet, please stop operation.
- When it is carried from cold environment to warm or humid environment, please do not use it immediately.
- DO NOT operate button on front panel with sharp things.
- High temperature or high pressure steam disinfection to the device is not permitted. Refer to User Manual in the relative chapter (7.1) for cleaning and disinfection.
- Do not have SPIROMETER immersed in liquid. When it needs cleaning, please wipe its surface with medical alcohol by soft material. Do not spray any liquid on the device directly.
- When cleaning the device with water, the temperature should be lower than 60°C.
- The display period of data is less than 5 seconds, which is changeable according to the end rate.
- When data can't be displayed at all times or other cases happened during testing, press "repeated measure" key to remeasure, or power off to restart.
- The device has normal life for three years since the first electrified use.
- When the data goes beyond the limits, the main screen shows "Error!".
- The device doesn't suit all users, if you can't get food measurement data, please stop using it.
- The device needs to be calibrated once per year or less.
- The device is forced SPIROMETER, according to the User Manual to use right to gain best result.

### EMC Declaration:

- When the device is installed or put into service, extra attention must be paid to EMC, as portable and mobile RF communications equipment with higher EMC interferences may affect this device.
- Internal components and cable should not be changed as this may decrease the IMMUNITY of the device.
- The SPIROMETER should not be used adjacent to or stacked with other equipment.

## Chapter 2. Overview

Forced Vital Capacity is the maximum expiration after taking a breath, it is an important test content in chest-lung diseases and respiratory health, and is indispensable in modern lung inspection projectors respiratory health, and is indispensable in modern pulmonary screening. At the same time, it has great significance in respiratory diseases, differential diagnosis, treatment evaluation and selection of surgical indications. Thus, with this rapid development of clinical respiratory psychology, clinical applications of lung capacity inspection are also gaining popularity.

Forced Vital Capacity is the maximum expiration after taking a breath, it is an important test content in chest-lung diseases and respiratory health, and is indispensable in modern lung inspection projectors respiratory health, and is indispensable in modern pulmonary screening. At the same time, it has great significance in respiratory diseases, differential diagnosis, treatment evaluation and selection of surgical indications. Thus, with this rapid development of clinical respiratory psychology, clinical applications of lung capacity inspection are also gaining popularity.

## 2.1 Characteristics

- (1). Ultra slim, concise and fashionable design.
- (2). Small in volume, light in weight and convenient to carry.
- (3). Low power consumption.
- (4). TFT display
- (5). Reflects lung function by measuring FVC, FEV1, PEF, etc.

## 2.2 Major applications and outreach

The SPIROMETER is a hand-held equipment for lung function testing. The product is suitable for hospital, clinic, family or ordinary test. The user is required to operate according to the user manual, no specialised training is needed, therefore the operation of the device is as simple and easy as possible.

## 2.3 Environment requirements

### Storage Environment:

Temperature: -40°C ~+60°C

Relative Humidity: 5%~95%

Atmospheric pressure range: 500hPa~1060hPa

### Workinh Environment:

Temperature: -40°C ~+60°C

Relative Humidity: ≤80%

Atmospheric pressure range : 700hPa~1060hPa

# Chapter 3. Principle

First, the examinee breathes in deeply, then puts his lips around the mouthpiece and blows all the air out as hard as he can, the exhaled gas is transformed into a rotating turbine airflow, then rotates the blade.

The exhaled gas is transformed into rotating air flow by turbine, then rotates the blade. The receiving part of the pair of infrared diodes (one for infrared emission, the other for reception) towards the blade is used to receive the infrared rays when the blade rotates. The received beam strength of the receiving diode will be different as the difference of blade angle, to form the varied signal in the same proportion in the receiving diode, which forms acquisition signal by SCM after processing. Finally, various measured parameters formed with the information processed by the microprocessor and displayed on the screen.

# Chapter 4. Technical specifications

## 4.1 Main performance

- Forced Vital Capacity (FVC).
- Forced Expiratory Volume in one second (FEV1) and judge the condition of the examinee by the ratio of FEV1 and FVC.
- Peak Expiratory Flow (PEF), 25% FVC flow (FEF25), 75% FVC flow (FEF75) and mean flow between 25% and 75% FVC (FEF2575).
- Display of flow-volume rate graph, volume-time graph.
- TFT display.
- Battery power display.
- Automatic shutdown when there is no operation within one minute.
- Has the functions of data memory, delete, load and review.
- Trend curve display.
- Scaling (Calibration)
- Rechargeable lithium battery.

## 4.2 Main Parameters

Maximum volume: 10L

Flow range: 2 L/s ~16 L/s

Volume accuracy: ±3% or 50mL (whichever is greater)

Flow accuracy: ±5% or 200 mL/s (whichever is greater)

### Classification:

EMC: Group I Class B.

According to MDD 93/42, classification of this medical device: IIa.


Type of protection against electric shock: Internally powered equipment.

Degree of protection against electric shock: Type BF applied part.

Degree of protection against water ingress: IPX0.

## Chapter 5. Installation

### 5.1 Front panel view

### 5.2 Installation

- (1). Take the turbine, align the arrowhead on the turbine with the triangle on the casing, inserted in the lower part, rotate it anticlockwise to close it. clockwise to close it.
- (2). Insert the disposable nozzle into the turbine inlet.

### 5.3 Accessories

- (1). User's Manual
- (2). USB data line
- (3). Disposable Mouthpiece
- (4). Power adapter (optional)
- (5). CD (PC software)
- (6). Nose clip (optional)

**Other types of adapters must meet the following conditions: output voltage: DC 5V; current output >500Ma, the power adapter must meet the requirements of the related EN60601 standards and be CE marked.**

## Chapter 6. Operation Guide



### 6.1 Method of application

#### 6.1.1 Installation

Take the turbine, align the arrowhead of the turbine with the triangle on the casing, inserted at the bottom, rotate it counterclockwise to close it, then insert the disposable nozzle into the inlet of the turbine.


#### 6.1.2 Measurement

- (1). Press and hold the "on" key to turn on the power after installation.
- (2). The device is in selection interface after power on as in Fig.2, press "up", "down" key to adjust, select "NO" for "Check" interface as in Fig. 3.
- (3). Then take maximum air and put your lips around the mouthpiece and blow as hard as you can, wait for a few seconds and the data will be obtained, the measurement is finished. (Note: "Yes" indicates that you can edit the patient information, exit after editing or return to "Checking" interface, the operation detail is as follows).

### 6.1.3 Parameter interface

In the test interface take the maximum air and put your lips around the nozzle and blow as hard as you can, wait a few seconds and the measurement parameter will come out automatically.

(Note: Status indicator bar: indicates the measurement status, displays the condition of the examinee by the ratio of the measured value and the predicted value. I.e. Compares the measured value with the reference value in the same situation, when the value is less than 50% indicates that it should be perceived and hospitalised in time; The value in the range 50%-80% indicates that it should be perceived: it is green when the value is higher than 80%, which is normal. The status value is optional, press "Control Settings" in the main menu and press "Denote Value" to select).  
 Flow rate-volume graph, volume-time graph: press "up", "down" key, two graphs will appear: Flow rate-volume and volume-time (as Fig.5-6). The three interfaces above are the main interface. In its interface press "menu" for main menu as Fig.7.


Fig.5 gráfico tasa de flujo-volumen


Fig.6 gráfico volumen-tiempo

### 6.1.4 Menu operations

By browsing "menu" in its main menu as Fig.7, the user can search for other parameters and can operate control settings, patient information, real time configuration, shutdown, etc. The detail of methodology is as follows:


Fig.7

(1). Other parameters

Press "Other Par" in the main menu in its submenu as Fig.8 which displays other parameters except the three parameters of the main menu (see the following for details), press the confirmation key to return to the main menu (pressing the "up", "down" key on its interface has no effect).

The test value of other parameters	
FEV1%	100%
FEF25	5.55
FEF 75	2.81
FEF2575	4.19

Fig.8

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
### a) Trend curve

Select “Trend Curve” to the trend curve determinant interface as Fig.10. Press “up” or “down” key to select determinant value, then press confirmation key to display the variation curve, convenient for user to contrast the data. If the data are too many, press “up” or “down” key on the curve interface to search all the trend curve data orderly. Press the confirmation key to exit the current interface and return to the control setting interface (as Fig.9).

Control Setting	
Trend Curve	
Review Fun	OPEN
Denote Value	
Delete Data	
Blue Tooth	CLOSE
Scaler	
Exit	

Please select determinant of trend curve	
FVC	
FEV1	
PEF	
FEV25	
FEV75	
FEV2575	

## b) Review Information


Fig.12

Select "Review Number" in its submenu. If the previous display status is "CLOSED" press the confirmation key to open the function (note: the review function can only be opened when the number of cases is more than one), the interface will jump to the Number interface as shown in Fig.12. In this interface use the "up" or "down" key to select the case number, press the confirmation key to review the selected case information. If the previous display is "OPEN" press the confirmation key to jump to the previous status (meaning: If the measurement has been completed before review, close this function, display "CLOSED" status, or jump to the measurement interface as shown in Fig.2).

**Attention:** When opening the review function, you can search all the measured data, the way is: Enter the main menu as Fig.7, select "Patient Info", press the confirmation key in its submenu as Fig.21, select "Number", press the confirmation key to go to the case number interface as Fig.12 (the function can only be selected in case of opening the review function), then use the "up" or "down" key to select the case number, press the confirmation key to search the corresponding information. If you want to search for other cases, repeat the above steps (If the case data is wrong, "Error!" will be displayed after pressing the confirmation key).

## c) Denoted Value

Select "Denoted Value" in its submenu as Fig.14, press the confirmation key to choose what will be the denoted value (as Fig.14 the denoted value is decided by the ratio of FVC and the mean value), after selecting press the confirmation key to exit the interface (Note: the status value is a percentage decided by the ratio of the measured value and the predicted value).


Fig.13


Fig.14

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
Fig. 20

f) Scaling Operation (Calibration)

Select "Scaler (Calibration)" in its submenu as Fig.17, the volume of Scaler (Calibration) can be selected, after selecting the volume it will enter the scaler (Calibration) interface as Fig.18. In this interface by moving the scaler once the interface will display "REPEAT", move the scaler continuously. After two successful times the scaler (Calibration) is ready, the interface will display "OK!". Finally the interface will jump to the previous interface before scaling (Calibration). If Fig.19 appears indicating that the scaler (Calibration) has error, repeat the scaling (Calibration) and until it is successful. If Fig.20 appears confirm the volume between the scaler (Calibration) and select accordingly until the scaling (Calibration) is completed successfully. If you need to stop the scaling (Calibration), press the confirmation key to exit to the previous interface before scaling (Calibration). (The interface before scaling (Calibration) refers to: If scaling (Calibration) has been scaled (Calibration) after the complete measurement will exit to the control setting interface, otherwise it will jump to the measurement interface as Fig.2).

g) exit

Select "Exit" to leave the "Control Configuration" interface and return to the main menu.

(3). Patient Information

Press "Patient Info" in its submenu as Fig.21 (Note: Fig.2 is the selection interface, and "Yes" indicates that you can edit the patient information).


Fig.21

a) Number

"Number" is the current patient data. Attention: If the review function is open, press the confirmation key at the case number interface to review the patient information, see "Review Information" in Section 6.1.4 for detailed operations.

b) Gender settings

Press "up" key, "down" key for "Gender", press key to select "FEMALE" OR "MALE".

c) Patient information setting (age, height, weight and nation)

Press "up", "down" key in the information interface for "Age", press "menu" for "set age" interface, press "up", "down" key to set age, press "up" or "down" key once increases or decreases a number to the ideal value. Press "menu" to interface "Patient Info". The setting operations "Height", "Weight", "Nation" are the same as "Age". The range of "Age" is 6-100, "Height" is 60-240cm, "Weight" is 15-250kg, select "Nation" according to the predicted standards, including "ERS, KNUDSON, USA".


Fig.22


Fig.23

d) Option Smoker, Drugs

Press “up”, “down” key in the patient information interface for “Nation”, press “down” key and “Smoker” appears, if you keep pressing “down” key, “Drugs” appears. Press the confirmation key, you can select “Yes” or “No” to edit “Smoker” or “Drugs”.

e) Exit

Press the “up”, “down” key to exit the “Patient Info” interface and return to the main menu.

#### (4). Switch off

Press “Power off” to turn off the device.


#### (5). Exit the main menu

In the main menu interface, press “Exit” to exit the main menu interface.

### 6.1.5 Repeated measurement

Press and hold the “repeat measurement” button for no more than 3 seconds to start another test, the information as in Fig. 12 will appear when the memory is full.

Press “up”, “down” key to set up, select “Yes” for “clear data” interface, when select, after clear data, exit to main interface to continue measurement, “NO” for main menu.

### 6.1.6 Loading

There are two types of loading methods:

- (1). Connect the device to the computer with data line, then the device will be in charging state.
- (2). Connect the device to the power supply with power adapter, then the device will be in charging state.

**DO NOT use the machine when it is in charging state with this power adaptor.**

**The red and green indicators are all highlighted in the charging state, the red indicator goes out when charging is complete.**

### 6.1.7 Load data

Install the PC software on the computer, the following figure will appear on completion.

Fig. 28

- (1). Connect the device with the computer by data line, double press the icon to open the PC software procedure.
- (2). Press the corresponding key to get upload data, delete cases, print information, history, select language, change PDF format, set examinee information, etc.
- (3). Press "exit" to exit the software, unplug the data line from the computer to finish uploading.

## 6.2 Attention

- Please check the device before use and confirm that it is working normally.  
it is working normally.
- Rechargeable lithium battery.
- It is recommended to measure the device in room.
- Excessive ambient light may affect measurement accuracy.  
Includes fluorescent lamps, dual ruby lights, infrared heaters, direct sunlight, etc.
- Intense subject activity or extreme electrosurgical interference may also affect accuracy.
- Please clean and disinfect the device after use according to the User Manual (7.1).

## Chapter 7. Maintenance, transport and storage

### 7.1 Cleaning and disinfection

Use medical alcohol to disinfect the device, let it dry or wipe it with a clean, soft cloth. It is necessary to clean the turbine periodically for accuracy, keep it clear of the luminous part and keep it away from dryers. Immerse the turbine in a detergent solution after use, clean it with clean water and let it dry vertically for a few minutes (but do not rinse directly with water), this type does not pollute the environment.

### 7.2 Maintenance

- (1). Please clean and disinfect the device before use according to the User Manual (7.1).
- (2). Please recharge the battery when the display shows low power (battery power is  ).
- (3). Recharge the battery soon after over-discharge. The device should be charged every six months if not used regularly. This guide can extend the life of the battery.
- (4). Please clean and disinfect the device after use to avoid infection. The device needs to be calibrated once a year (or according to the hospital calibration schedule).

It can also be performed at the state designated agent or simply contact us for calibration.

### 7.3 Transport and storage

- (1). The packaged device may be transported by ordinary transport or in accordance with the contract of carriage. The device may not be transported mixed with toxic, harmful, corrosive material.
- (2). The packaged device must be stored in a room without corrosive gases and with good ventilation. Temperature: -40°C~60°C; Humidity: <95%.

## Chapter 8. Troubleshooting

Problem	Possible reason	Solution
The device cannot finish the measurement for a long time and the data is not displayed.	Start-up speed is too low, device does not measure.	Re-measure according to the user manual..
	Device malfunction	Press the "repeat measurement" key to re-measure or switch off to restart.

The figures are incorrect and out of order.	The power has been switched off normally.	Delete the current case and re-measure
	The operation is wrong	Operate normally according to the user's manual.
	Device malfunction	Please contact your local service centre.
The device cannot be switched on	The battery is empty or almost empty.	Please charge the battery.
	Device malfunction	Please contact your local service centre.
The visualisation suddenly disappears.	The device shuts down automatically when there is no operation for one minute.	Normal
	The battery is empty or almost empty.	Please charge the battery
The device cannot be used full time after charging.	The battery is not fully charged.	Please charge the battery
	The battery is broken	Please contact your local service centre.
The device cannot be used full time after charging.	The battery is broken	Please contact your local service centre.

## Chapet 9. Symbol Key

Symbol	Meaning
	Caution, see instructions for use.
	WEEE (2002/96/EC)

	Type BF Applied part
	Full energy
	Low energy
Error	Measured value is beyond limits
	DO NOT REUSE
	Status indicator bar

## Chapter 10. Entering parameters

### Measured parameters

Parameter	Description	Unit
<b>FVC</b>	Forced Vital Capacity	L
<b>FEV1</b>	Forced Expiration Volume in one second	L
<b>PEF</b>	Peak expiratory flow	L/min
<b>FEV1%</b>	$FEV1/FVC \times 100$	%
<b>FEF25</b>	25% FVC flow	L/s

<b>FVC</b>	Mean flow between 25% and 75% of FVC	L/s
<b>FEV1</b>	75% FVC flow	L/s

CONTEC MEDICAL SYSTEMS CO., LTD  
 No. 112 Qinhuang West Street, Economic &  
 Technical Development Zone, Qinhuangdao,  
 Hebei Province, P.R.C.  
[cms@contecmed.com.cn](mailto:cms@contecmed.com.cn)

PROLINX GMBH  
 Brehmstr, 56, 40239, Dusseldorf Germany  
[med@eulinx.eu](mailto:med@eulinx.eu)

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GRUPO R. QUERALTÓ S.A.