MANUAL DE INSTRUCCIONES

ELECTROCARDIÓGRAFO


ECG600G

Manual de instrucciones

ES

Instruction manual

EN

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
ECG600G

NOTICE

Notice to the user and/or patient; any serious incident that has occurred in relation to the product must be informed to the manufacturer and to the local authority of the Member State in which the user and/or patient is established.

INDICATIONS

- In any type of patient to detect episodes of hypoxaemia difficult to detect clinically.
- Situations of respiratory failure of any aetiology.
- Emergency tracheal intubation and connection to mechanical ventilation.
- As a heart rate monitor in any emergency procedure.

1. INTRODUCTION

This product is an electrocardiograph with 12 simultaneous ECG leads that allows recording ECG waveforms with a thermal printing system. Its functions are as follows: recording and displaying ECG waveforms in automatic manual mode; self-measuring and self-diagnosis of ECG waveform parameters; indication of lead off and paper out; language switching (Chinese/English); AC or DC power supply; rhythm lead can be selected to detect abnormal heart rhythms; database management, High resolution 800*600 LCD display, 1728 dots wide printing, can print twelve-channel ECG waveforms. Operated by keyboard or touch screen.

1.2. Technical specifications

- Environmental conditions
 - Operation
 - a) Environment temperature 5°C~40°C
 - b) Relative humidity: 25%~95% (no condensation).
 - c) Power supply : AC: 100V~240V, 50/60Hz. Rechargeable lithium battery 7.4V, 3500 mAh
 - d) Atmospheric pressure: 700 hPa~1060 hPa
 - Transport and storage:
 - a) Ambient temperature: -20°C~55°C
 - b) Relative humidity: ≤95%.
 - c) Atmospheric pressure: 500 hPa~1060 hPa
- Input mode: Floating circuit and defibrillation protection
- Lead: 12 standard leads,
- Patient leakage current: <10μA
- Input impedance: ≥2.5 MΩ
- Frequency response:

Rated input amplitude	Input frequency and waveform	Relative output response
1.0	0.67 Hz~40 Hz, Sine wave	±10% ^a
0.5	40 Hz~100 Hz, Sine wave	+10%, -30% ^a
0.25	100 Hz~150 Hz, Sine wave	+10%, -30% ^a
0.5	150 Hz~500 Hz, Sine wave	+10%, -100% ^a
1.5	≤1 Hz, 200 ms, Triangle wave	+0%, -10% ^b
^a Relative to 10 Hz ^b Relative to 200 ms		

- Time constant: >3.2s
- CMRR: >105 dB
- Filter: power frequency (AC 50/60 Hz), myoelectricity (15 Hz/35 Hz (-3 dB)), baseline drift filter.
- Recording mode: thermal printing system
- Printing paper specifications: 110 mm (width)*20 m (height), high speed thermal paper.
- Paper speed: 12.5 mm/s, 25 mm/s, 50 mm/s, error ±5%
- Sensitivity options: 5, 10, 20 mm/mV, accuracy is ±2%; Standard sensitivity is 10 mm/mV ±0.2 mm/mV.
- Automatic recording: works according to recording format and automatic mode, automatic lead change, automatic measurements and analysis.
- Pace recording: operates according to pacing format and pacing mode, automatic measurements and analysis.
- Manual recording: works according to recording format, manual lead change.
- Measurement parameters: HR, P-R interval, P duration, QRS duration, QRS duration, T duration, Q-T interleave, Q-Tc, QRS axis, T axis, R (VS), S (V1), R (VS)+S(V1)

- Product safety type: Class I type CF defibrillation-proof applied part
- Polarization resistance voltage: ± 610 mV
- Noise level: ≤ 12 μ Vp-p.
- ECG signal input sampling frequency: 32 kHz
- Waveform data processing sampling frequency: 1 kHz
- Sampling precision: 24-bit
- Minimum detection signal: 10 Hz, 20 μ V (peak-peak value) deflected sinusoidal signal can be detected.
- Pacing detection channel: standard II
- Accuracy of input signal: $\pm 5\%$
- Amplitude quantization: $\leq 5\mu$ V/LSB
- Dimensions: 315 mm (length)*215 mm (width)*92 mm (height).
- Net weight: 1.7 kg
- Interchannel time deviation: <100 μ s

1.3 Main Characteristics

- 7-inch high-resolution color LCD display, operated either through the touch screen or the function buttons.
- Sync collection for 12-lead ECG, uses digital signal processing technology to conduct AC filter, baseline filter and EMG filter on ECG signals, in order to get high-quality ECGs.
- Display of 3/6/12-lead ECG on one screen, and print mode, sensitivity, paper feed, filter state and other information, which facilitates comparative diagnosis.
- The device can be powered either by AC or DC (can adapt to 50/60 Hz AC frequency), with built-in rechargeable lithium battery and charging circuit, perfect battery overcurrent and overvoltage protection circuit.
- Multiple print mode and format, including automatic 12x1 6x2+1 (rhythm lead), 6x2, 3x4+2 (rhythm lead), rhythm 12, rhythm 10, rhythm 8, rhythm 6 and manual Trend graph and histogram of RR interval can be printed. Printed waveform length is adjustable. With timed print function to satisfy various application requirements.
- Rhythm leads can be arbitrarily selected to facilitate the observation of abnormal heart rate.
- Clinical information such as patient's number, name, gender, age, height and weight can be input.
- The built-in large capacity memory can store 1000 medical records, making it easy for the physician to review medical records and statistical information.
- Multi-language interface and report.

1.4 Software Overview

The ECG analysis program shows the results after analyzing the form of the electrocardiogram, providing auxiliary reference for doctors to make a diagnosis. The analysis results cannot be used as the only standard for a diagnosis. A comprehensive evaluation should be made by professional electrocardiogram technicians and physicians according to clinical experience and other test results.

The device is intended for use on all patient populations, which is decided by the clinical doctor. The analysis program only provides ECG analysis for patients above 3 years old (included).

- Name of software: ECG600G embedded software.
- Software specification: none
- Software version: V1.6.7
- Version naming rules: V <major version number>.<minor version number>.<revision version number>
- The version of the software can be obtained in "About".
- Name: ECG algorithm.
- Type: Mature algorithm.
- Use: to convert ECG signals of human body into intuitive waveform images and the analyzing.
- Clinical function: An electrocardiogram is an important method for clinical diagnosis of cardiovascular disease. How to use computer quickly, automatically and accurately analyze ECG has been a relevant topic among scholars. The ECG algorithm is key to determine the effectiveness of diagnosis and treatment of patients with heart disease.

2. SAFETY PRECAUTIONS

- Ensure that the device is placed on a flat level work table. Avoid strong vibrations or impacts when moving it.
- When working with AC power, the power cord must be 3-core, the frequency and voltage value of the AC power source must match the identification on the manual and have sufficient capacity. When the provided three-core power cord cannot be used, please use the built-in DC power supply or replace the three-core power cord that meets the standard requirements.
- A perfect power supply system and grounding are necessary in the room.

Warning: to avoid the risk of electric shock, the device must be connected to a power supply with protective grounding.

- If you have any doubts regarding the reliability of the protective grounding cable, or the connection can not be guaranteed, the

device must be run with built-in DC power supply.

- The safety requirements have been fully considered in product designing, but the operator can't ignore the observation of the patient and device. Cut off the power or take off the electrode when necessary to ensure the patient's safety.
- Please turn off the device and unplug the power cord before replacing the fuse or cleaning and disinfecting. Don't use sharp materials to rub the screen.
- Keep the device away from water, don't use or store it in places with high air pressure, humidity or temperatures over the recommended conditions, bad ventilation or too much dust.
- Do not use the device in places with flammable anesthetic gases or other flammable chemicals as there is a danger of explosion or fire.
- Do not use the device in medical hyperbaric oxygen chambers, as it may cause an explosion or start a fire.
- The device is not intended to act directly on the human heart. If this device is used with a cardiac defibrillator or other electric stimulating devices at the same time, single-use electrodes and ECG lead cables with defibrillation-proof function should be selected. It is better not to use this device with other electric stimulating devices at the same time. If it is necessary, there must be a professional technician present to provide guidance and the selected accessories should be approved by the manufacturer.

Warning: Do not operate the instrument on parts of the human body with wounds and do not perform measurements on parts with surface wounds.

- When the electrocardiograph is used together with a high-frequency electrosurgical knife, the ECG electrode should be kept away from the knife to avoid burns or burning of the electrode wires due to high-frequency sparks.
- When the electrocardiograph is used together with a defibrillator, the operator should avoid contact with the patient or the bed. The defibrillation electrode should not touch the ECG electrode directly to avoid burns to the device or the patient.
- Please do not use the electrocardiograph in environments interfered by high-power devices such as high voltage cables, X-rays, or ultrasonic machines. Keep the device away from emission sources such as mobile phones.
- If other equipment is connected with the ECG device, it must be a class I device that complies with IEC 60601-1. Because the total leakage current may hurt the patient, the monitoring of leakage current is carried out and taken charge by the connected equipment.
- Notes related to EMC:

The device complies with the safety standards for medical electrical equipment or system electromagnetic equipment compatibility in IEC60601-1-2. Electromagnetic environments exceeding the IEC60601-1-2 standard may cause harmful interference to the device or prevent the device from performing its intended function or degrade its performance. Therefore, if there is a phenomenon that does not match its function during use, be sure to confirm and eliminate adverse effects before continuing use. Corresponding precautions for this situation are given in this manual.

The device or system should not be used near or stacked with other devices. If it must be used near or stacked with other devices, it should be observed and verified that the device is working normally under the configuration used.

Use of accessories, transducers and cables other than those specified by the manufacturer of the device or system as spare parts for the internal components may result in increased emissions of the devices or system and reduced immunity,

- Effects from radiated electromagnetic waves:

The use of a mobile phone may affect the operation of the device. When installing medical electrical equipment, be sure to remind people around the device to turn off mobile phones and small radios.

- Effects from shock and conduction electromagnetic waves:

High frequency noise from other equipment can enter the device through the AC socket. Please identify the source of noise, if possible, stop using the equipment. If the equipment can not be deactivated, use noise cancellation equipment or take other measures to reduce the impact.

- Effect from static electricity:

Static electricity in a dry environment (indoors) may affect the operation of the device, especially in winter. Before using the device, humidify the indoor air or discharge the static electricity from the cable and operator.

- Effect from thunder and lightning:

If there is lightning nearby, it may cause a voltage surge in the device. If you are concerned about danger, disconnect the AC power and use the internal power supply.

• Notes concerning ECG waveform measurement and analysis:

- The identification of P wave and Q wave is not always reliable with intensive EMG or C interference. Neither are the ST segment and T wave with baseline drift.

- Winding and unclear end position of S wave and T wave may cause error in measurement.

- When R wave is uninspected caused by some leads off or QRS wave low voltage, the heart rate measurement may deviate greatly from the correct one.

- In case of QRS low voltage, ECG axis calculation and border-point identify of QRS wave are not always reliable.

- Occasionally, frequent ventricular premature complexes may be identified as dominant beat.

- Merging of versatile arrhythmia may result in unreliable measurements because of the difficulty in distinguishing P wave in such situation.

- The device has an automatic analysis function that analyzes the obtained ECG waveform without reflecting all the patient's status. The results of the analysis may sometimes not comply with the physician's diagnosis. Therefore, the final conclusion needs to be comprehensively analyzed by doctor in combination with analysis results, patient clinical characterization and other test results.

3. WORKING PRINCIPLE AND STRUCTURAL CHARACTERISTICS

3.1 Working principle and its block diagram

3.1.1 The power supply unit

- Principle of power supply

After the AC power supply enters the switching power supply, it is converted to 12V DC voltage and supplied to the DC-DC power board, it also provides constant voltage current limiting charging for the rechargeable lithium battery in the device through the DC-DC circuit and generates +5V and +3.3V voltage through the power conversion to supply power to the corresponding modules. At the same time, the lithium battery in the device can independently satisfy working requirements of each module in the device through the buck-boost circuit.

WARNING: The principle block diagram and component list are only available to service stations or maintenance personnel designated by our company.

3.1.2. Signal acquisition unit

The signal acquisition unit uses a floating setting, which is a signal acquisition and processing system, including analog circuit part and A/D conversion (with sampling accuracy of 14 bits) and data processing part. The analog circuit consists of signal following, simplification, anti-aliasing low-pass filtering, lead-off detection and overload detection. CPU system is responsible for coordinating the work of each circuit such as the A/D converter, the lead-off detection circuit and the overload detection circuit, in order to achieve signal acquisitions, processing and lead-off detection. Control information and A/D conversion and data acquisition between the floating circuit and the solid circuit are transmitted through the optoelectronic coupler.

3.1.3 Control unit

- Principle of control unit

The control system consists of printing system, button system, liquid crystal display system and signal acquisition system. The ECG signal sent from the signal acquisition system through the high-speed optoelectronic couple is received by the CPU system, after digital filtering, gain adjustment and motor drive, it is sent to the printing system to print the ECG waveform. After the printing is completed, the CPU system processes waveform measurement and analysis. The CPU system also receives an interrupt signal and button code from the button system to complete the interrupt processing. In addition, the lead-off signal, paper out detection, battery voltage management and automatic power-off are also managed by the CPU system. The liquid crystal controller receives data and commands from the CPU system to complete the display of the control state of the device.

- Principle block diagram is shown in Figure 4-1


Figure 4-1 Block diagram of control unit

3.2 Name of parts and functions

3.2.1. Front view


Figure 4-2 Front view

1. Paper compartment cover. Keeps the paper compartment closed, holds the printing paper.
2. Display screen. Displays patient's ECG and related information.
3. Button area. Controls the operations of the device and enter the information.

WARNING:

- Do not put heavy objects on the screen or hit it, otherwise the screen will be damaged.
- If the device is not in use, cover it to prevent spills on the screen.
- Do not use sharp objects to operate the buttons, otherwise it may cause permanent damage to the buttons.

3.2.2 Side view


1. Input socket: connect with AC power cord.
2. Equipotential terminal: Connect with the potential equalization conductor.
3. Lead cable interface: Connect with lead cables.
4. USB interface: Communicate with the computer. The ECG data can be transmitted to a computer. Many functions can be achieved with its use such as archiving, managing and analyzing ECG data which facilitates clinical research organization teaching and training.

WARNING:

- Lead cables must be disconnected from patient before connecting with computer via the USB interface.
- Operator must not touch the USB interface and patient at the same time.



3.2.3. Buttons

1. Startup indicator: lights in green after turning on the device.
2. Power status indicator: Green indicates that the AC power supply is being used. At this time there is no battery in the device or the battery is full. Red and green alternating means the battery is being charged.
3. Direction buttons: including up, down, left, right and OK buttons, quick and convenient.
4. PRINT: Print the sampled ECG waveform or end the printing.
5. MODE: When the device is in sampling interface, use the MODE button to select the print mode.
6. Setting button for acquisition system: Collect ECG waveform and set the display mode.
7. ON/OFF: When the device is turned on, short press this button and it will prompt whether to shut down the device. Long press this button to turn off the device.

3.2.4 Symbols

~AC	Alternating current
	Equipotential point, the equipotential point of this device is combined with the protective grounding.
	Caution! Refer to the accompanying document.
	Type CF equipment, with defibrillation-proof function
	USB interface
◀ PATIENT	Lead cable socket
	Serial number
	Manufacturer
	Date of manufacture
	Batch code
	Latex free
	Atmospheric pressure limitation

	Temperature limitation
	Humidity limitation
	This way up
	Fragile, handle with care
	Keep away from rain
	Stacking limit by number
	This item is compliant with Medical Device Directive 93/42EEC of June 14, 1993, a directive of the European Economic Community
	Refer to instruction manual/booklet
	Catalogue number
	Authorized representative in the European Community

4. OPERATION PRECAUTIONS

4.1 Precautions before use

- For safe and effective use, please read the user manual carefully before operation.
- Check to ensure that the device is in good condition.
- The device should be placed on a flat surface and moved gently to avoid strong vibrations or shock.
- Check to ensure that the lead cables are correctly connected and the device grounding is correct.
- The AC frequency and voltage should comply with the requirements and enough current capacity should be guaranteed.
- When using the battery for power supply, check to ensure that the battery voltage and battery status are in good condition and the battery has enough power.
- When the device is used together with other equipment, all devices and equipment should have equipotential grounding in order to protect the user and the operator.
- Install the device where it can easily be grounded in the room. Do not allow the patient and patient-connected lead cables to come into contact with other conductor parts, including the earth or a hospital bed.
- Clean the lead cable with neutral solvent. Do not use alcohol-based cleaners or germicides.
- Ensure that the device is running within the normal environmental temperature range of 5°C-40°C. If the device is stored at a higher or lower temperature, leave it in the operating environment for approximately 10 minutes before use in order to ensure correct operation.

4.2 Precautions during use

- The printing can be started after the ECG waveform is stable.
- During use, the doctor should observe the patient carefully and should not leave the operating site. If necessary, turn off the power or remove the electrode to ensure patient safety.
- The patient and the device can only be connected via lead cables through the electrodes, in order to avoid the patient touching other parts of the device or conductor.
- The patient should not move during operation.
- Maintenance and repair to the device or accessories is not allowed during use.

4.3 Precautions after use

- Set the states of all functions to initial states.
- Cut off the power, gently remove the electrodes and limb clips, then remove the lead cables, do not pull with force.
- Clean the device and all accessories and store them for the next use.

5. PREPARTIONS BEFORE OPERATION

5.1 Installation of recording paper

- The device adopts high-speed recording paper, its specification is 110 mm(W) x 20 m (L).
- The installation method of recording paper is described as below:
 - Slide the cover switch to the left to open the paper compartment cover. Take out the paper axis, insert it into the paper roll. The paper side with grids should be facing downwards. Then install it to the correct position in the paper compartment.

Figure 5-1 Installation of recording paper

- Pull out the recording paper from the slot of the paper compartment cover and close it.
- If the recording paper runs out during recording, the device will stop printing automatically and the screen will display a prompt of lack of paper, as shown in figure 6-2.

Figure 5-2 Lack of paper prompt

5.2 Power supply connection

5.2.1 AC

Insert one end of the provided three-core power cord into the device's input socket and insert the other end into a three-core power socket that meets the requirements. Ensure that the connection is secure and reliable and the device is automatically grounded. When the device is used in conjunction with other medical equipment, use the supplied potential equalization wire to connect the equipotential terminal of the device to the equipotential terminal of the connected equipment to prevent leakage current and protect the device.

5.2.2 Battery

The device has a built-in rechargeable lithium battery, which does not need to be re-installed by the user. Check the battery's power and status before use.

Warning: Connect one end of the potential equalization wire to the equipotential terminal of the device and connect the other end to the ground to enhance the reliability of the grounding. Do not use other pipes as ground wire, otherwise, the patient may be in danger of electric shock.

5.3. Lead cable connection

Connect the lead cable to the lead cable interface on the device and fasten it to the device with the fixing knobs at both sides of lead cable in order to prevent bad connection and affecting the detection.

Warning: The lead cable interface can not be used for other purposes except as the input interface of ECG signals.

5.4 Electrode installation

Proper installation of the electrodes is an important part of accurately recording the electrocardiogram. Make sure the electrodes are in good contact. Old and new electrodes or reusable electrodes and disposable electrodes cannot be used at the same time. If different types of electrodes are used together, it may cause a longer polarization time, which will seriously affect the ECG recording. The electrode or lead plug must not touch other object surfaces or conductors, such as metal beds. Please replace them all when updating the electrodes.

Warning: do not test on wounded skin.

5.4.1. Chest electrodes

As show in Figure 6-3:

The chest electrodes should be installed to the following parts:

C1 (V1): The fourth intercostal space at the right sternal margin.

C2 (V2): The fourth intercostal space at the left sternal margin.

C3 (V3): Between C2 and C4.

C4 (V4): The intersection between the midclavicular line and the fifth intercostal space.

C5 (V5): Left anterior axillary line on the same plane as C4.

C6 (V6): Left midaxillary line on the same place as C4.

Clean the chest skin where the electrodes will be installed with alcohol and apply some conductive paste to the skin (about 25 mm in diameter) and the edge of the chest electrode suction cup. Squeeze the suction ball to install the chest electrodes in positions C1-C6.


Figure 5-3 Installation of chest electrode

Warning: the conductive paste coating should be separated from each other and the chest electrodes should not touch each other to avoid short circuit.

Warning: please use qualified conductive past to avoid damaging the skin.

5.4.2 Limb electrodes

The limb electrodes should be placed on the skin of both hands and feet. Before connecting, clean the skin of the electrode installation area with alcohol, then apply a small amount of conductive paste on the clean skin. The electrode connection of the limbs is shown in figure 5-4.


Figure 5-4 Installation of limb electrodes


5.4.3. Colors of lead cables

Electrode position	European standard		American standard	
	Mark	Color	Mark	Color
Right arm	R	Red	RA	White
Left arm	L	Yellow	LA	Black
Left leg	F	Green	LL	Red
Right leg	N/RF	Black	RL	Green
Chest 1	C1	Red	V1	Red
Chest 2	C2	Yellow	V2	Yellow
Chest 3	C3	Green	V3	Green
Chest 4	C4	Brown	V4	Blue
Chest 5	C5	Black	V5	Orange
Chest 6	C6	Purple	V6	Purple

Warning

- It is recommended to install the lead cables after turning off the device.
- Apply appropriate amount of conductive paste on the electrode when installing the electrode.
- If the ECG waveform does not appear for a long time, check if the electrode is in contact with the skin.

5.4.4 Lead method and system

5.4.5 Lead-off and overload indication


The device can check the connection status of the lead at any time. If lead-off or overload is detected, the screen will display corresponding lead code on the top left corner, as shown in figure 6-2.

Warning:

- In the lead off prompt area, red font represents lead-off, yellow font represent overload.
- When the connection between lead cable and patient or the device is not reliable and the ECG signal can not correctly transmitted, the device displays lead-off.
- In the printed report, lead-off is marked with "°", and lead overload is marked with "+".

6. OPERATION INSTRUCTIONS AND PARAMETER SETTING

6.1 Main interface

A. Time: System time can be set in  so the detail time of the ECG recording can be noted.

B. Memory Usage: It directly shows the capacity of memory space according the actual usage. The green part represents the memory usage, the white part represents the remaining memory space.

C. Total number of cases stores in memory

D. Battery level (refer to 8.1)

E. Functional panel

Select to enter the sampling interface. Generally, the device will automatically enter this interface after powering on.

Select to enter the case management interface. In this interface the user can query, modify or delete case information.

Select to view the placement of leads.


Select to set time and date.

Select to perform system setting.

Select to perform sampling setting.

Select to set the parameters used in automatic analysis.


Select to set the print mode, print style and print content, etc.

Select to view the information about our company, software version.

Quick switch: use the     buttons on the keyboard to switch between each function module, after selecting press  button to enter its setting.

Quick setup: Click the functional module on the screen to quickly set the corresponding function.


6.2 Sampling interface

Click  on the main interface or press the  button to enter the sampling interface.

Warning: There is input time of cases in the system setting, therefore, in default, case information should be input before sampling.

The sampling interface provides several lead display modes, including 3-lead, 6-lead and 12-lead. The following figure uses 12-lead as an example:

- Status bar
 1. HR: current sampled heart rate value.
 2. Lead-off and overload: In demo mode, it displays "DEMO ECG". In sampling mode, it displays the detected lead status. A red lead icon represents lead-off. A yellow lead icon represents overload.
 3. Print status indication:

Display content	Explanation
Process...	It is printing.
Waiting...	It is finishing printing.
No Paper	Lack of paper, user should restart printing after loading paper.
Print timeout	Communication failure between this system and printing sub-system.
ECG timeout	Communications failure between this system and sampling sub-system.
Low power	Low power, it cannot start printing.

• Display field

The screen displays sampled 12-lead ECG waveform, by double clicking the waveform, you can switch between 3-lead, 6-lead and 12-lead.

• Operation field

Control the print display mode of the device through corresponding operation settings.

1. Speed: click the screen to switch the speed between 12.5 mm/s, 25 mm/s and 50 mm/s.

2. Gain (sensitivity): click the screen to switch the gain between 5 mm/mV, 10 mm/mV and 20 mm/mV. The calibration voltage signal and overall gain (sensitivity) can be checked by calibration function.

3. Filter: click the screen to switch the filter between no filter, AC filter, EMG filter, and DFT filter.

4. Print mode: use the  button to switch the print mode between Manual, Auto 6x2+1, Auto 6x2, Auto 4x3+1, Auto 4x3, Auto 3x4+1, Auto 3x4, Auto 2x6+1, Auto 2x6, Rhythm 7, Rhythm 6, Rhythm 5, Rhythm 4, Rhythm 3 and Rhythm 2.


5. Print/End print: use the  button to start or end the printing operation.

- Auto mode: After starting to print, the system automatically prints and stores the real-time 12-lead ECG waveform. The length is determined by the relevant settings in the print setup. Based on the settings, the automatic analysis data and conclusions are printed and the system automatically ends printing.
- Manual mode: After starting to print, the user needs to switch the lead to print the waveform of different leads, that is, the ECG printed in the manual mode is asynchronous and the data is not saved. The user needs to press the PRINT button again when the print needs to be terminated.


6. End sampling: After the device starts sampling, use the  button to end the sampling and back to the main interface.

6.3 Case information input interface


Due to the difference of system setup (refer to 6.6), users can choose to input case information before or after the sampling or not input the case information. The dialogue box is shown as below.

After selecting a box to edit, press the  button to show the keyboard as shown below. Selecting "Caps" can switch between numbers, lowercase letters, uppercase letters and symbols. "Space" is the space key, press it to enter a space; "Backspace" is the backspace key, press it to delete the last character entered. Click "OK" to confirm the entry and exit the interface.


The keyboard may have input restrictions according to the limitation of content. The restricted keys will be grayed out and unavailable, as shown below:

6.4 Case management

In the main interface, click  to enter the case management interface, as shown below:

The above interface shows all medical records stored in the device. User can search necessary cases by query function in the interface (refer to 7.4.1), modify or delete case information by edit function and review stored case information.

Click  to jump to the first page of case list.

Click  to jump to the last page of case list.

Click  to jump to the previous page

Click  to jump to the next page.

In the operation field, "List All" contains menus of export cases in the device to the file (cms60og_archive) of a U-disk via the USB interface.

6.4.1 Query


Click "Query" to enter the query interface shown below. Input query conditions and click "Query" to obtain expected results. After clicking "Clear", the system will clear all query conditions that have been entered.

"Cond.And" and "Cond.Or" indicate the matching mode of the query conditions. You can choose one of the two. If you select "Cond.And", the displayed query results will satisfy all the input conditions at the same time; if you select "Cond.Or", the displayed query results only need to meet any of the conditions entered.


Suggestion: When there are many cases, it would be better to input accurate query conditions and choose "Cond.And" to quickly find the case.

6.4.2 Review

In the case management interface, select a case to be reviewed, click "Review" to display the case information. The user is allowed to modify the patient's information, after clicking "Save", the information will be changed. Please note that the modification is irreversible. Make sure the input information is correct, click "Review" to enter the review interface, which is similar to the sampling interface.

Click "Report" to enter the diagnosis interface shown below, double-clicking to view a single QRS waveform. If a case length is larger than 6 seconds, clicking "Start" could review it in dynamic mode and the button turns to stop during dynamic reviewing. "<" and ">" and "Report buttons" are unavailable; if case length is shorter than 6s, "Start" is gray and unavailable.

In review interface, users can use  button to switch print mode and use  to print.

6.5 Date and time setup

In the main interface, click  to enter the following interface to set date and time.


In the current interface, the user can switch between items using the buttons   and adjust the content of the item with the buttons  . It can also be changed through the touch screen function which is convenient and quick.

6.6. System setup

In the main interface, click  to enter the system setup interface, as shown below:

After clicking "Default", the system will restore all settings to default.


The optional content of each setting item and its description are show in the following table:

Item	Options	Description
Back-light	[30 seconds]/[1 minute]/[2 minutes]/[5 minutes]/[10 minutes]/[Always on]	If there is no operation after reaching the set time, screen backlight will turn off. If it is set to "Always On", the backlight will always stay on.
Light-degree	[40% degree]/[60% degree]/[80% degree]/[100% degree]	After setting light degree, the screen will display different backlight strength.
Auto Off	[None]/[1 minute]/[3 minute]/[5 minutes]/[10 minutes]/[15 minutes]/[30 minutes]/[60 minutes]	If there is no operation after reaching the set time, the system will automatically turn off. If it is set to "None", the system will always keep on.
Low Power	[None]/[Only once]/[Always]	It determines the alarm method that the device uses with low power.
Filter Freq	[50 Hz/35 Hz]/[50Hz/25Hz]/[60Hz/25Hz]/[60Hz/35Hz]	To set the parameters of AC filter and EMG filter.
Language	[English]/[Chinese], etc.	To set the default language of system.
K-B Sound	On/Off	If it is selected, the button makes sound while pressing, otherwise, there will be no sound.
Demo Mode	On/Off	If it is selected, the system will run in Demo mode; otherwise, the system will run in sampling mode.

6.7 Sampling setup

In the main interface, click  to enter the sampling setup interface, as shown below:

After clicking "Default", the system will restore all settings to default. The optional content of each setting item and its description are shown in the following table:

Item	Options	Description
AC Filter	[ON]/[OFF]	Turn on or off the AC filter.
EMG Filter	[ON]/[OFF]	Turn on or off the EMG filter.
DFT Filter	[ON]/[OFF]	Turn on or off the Baseline filter.
Info Input	[Before]/[After]/[NoneInput]	Set to input the case information before or after sampling or no input
Show style	[3 leads]/[6 leads]/[12 leads]	Set the display method of ECG
Show gain	[5 mm/mV]/[10 mm/mV]/[20 mm/mV]	Set the gain of displayed ECG
Show speed	[12.5 mm/s]/[25 mm/s]/[50 mm/s]	Set the speed of displayed ECG, during printing, the auto mode and rhythm mode does not support the speed of 5 mm/s, 10 mm/s or 12.5 mm/s

6.8 Print Setup

In the main interface, click  to enter the print setup interface, as show below:

After clicking "Default", the system will restore all settings to default.


The optional content of each setting item and its description are shown in the following table:

Item	Options	Description
Print Mode	[Auto 6x2+1]/[Auto 6x2]/ [Auto 4x3+1]/[Auto 4x3]/ [Auto 3x4+1]/[Auto 3x4]/ [Auto 2x6+1]/[Auto 2x6]/ [Rhythm 7]/[Rhythm 6]/ [Rhythm 5]/[Rhythm 4]/ [Rhythm 3] and [Rhythm manual]	The system takes the selected option as default print mode.
Lead Gain	[Smart]/[Current]	The option selected will be used as printing gain mode. "Smart" means the system will adjust gain automatically to fit paper height; "Current" means it will use screen waveform gain as the of printing.
Auto Strip	[3 sec]/ [4 sec]/ [5 sec]/ [6 sec]/ [8 sec]/ [10 sec]/ [15 sec]/[20 sec]/[25 sec]	The system takes the selected option as the print time length of each strip.
Rhythm Strip	[10 sec]/ [15 sec]/[20 sec]/[25 sec]/[30 sec]	When "Print Mode" is set to "Rhythm 12" "Rhythm 10" "Rhythm 8" or "Rhythm 6", the system takes the selected option as the print time length of each waveform.
Average QRS	[6x2+Mark]/[6x2]/[4x3+Mark]/ [4x3]/[3x4+Mark]/[3x4]/[None]	When "Print Mode" is set to "Auto" or "Rhythm", the system uses the selected format to print the average QRS waveform
Auto-Diag	[All]/[Data]/[Conclusion]/ [None]	The diagnosis contains data and conclusion two parts, which can be chosen by user as needed.
Periodic	[OFF]/[per 1 min]/[per 2 min]/ [per 5 min]/[per 10 min]/[per 20 min]/[per 30 min]/[per 60 min]	During the ECG acquisition process, the system will automatically activate the printing operation according to the selected time interval. When the printing mode is manual mode, the printing will output "Auto 12x1" format, otherwise, it will output according to the current setting mode.


Note 1: The settings of auto strip, rhythm strip, average QRS, auto diagnosis and periodic print are only optional in auto mode and rhythm mode.

Note 2: If print time length is less than 8s, the sampling and analysis time will be 8s; if the print time length is equal to or greater than 8s, the sampling and analysis time keeps the same print time.


6.9 Lead placement

In the main interface, click  to view the schematic diagram of lead placement, the European Standard is shown as below:


The american standard is shown as below:

6.10 About

In the main interface, click  to view the information regarding the device, as shown below:


7. TROUBLESHOOTING

7.1 Auto shutdown

- The battery is almost running out, which causes overdischarge protection circuit action.
- The voltage of AC power supply is too high, which causes overvoltage protection circuit action.


7.2 AC Interface

- Is the device grounded reliably?
- Are the electrodes or lead cables connected correctly?
- Are the electrodes and skin daubed with enough conductive paste?
- Is the metal bed grounded reliably?
- Is the patient touching the wall or metal parts of the bed?
- Is the patient touching other people?
- Is there a high power electric equipment operating nearby such as an X-ray machine or ultrasonic device?

Warning: If the interference can not be removed after taking the above measures, please use an AC filter.


7.3 EMG Interference

- Is the room comfortable?
- Is the patient nervous?
- Is the bed too narrow?
- Is the patient speaking during the recording?
- Are the limb electrodes too tight?

Warning: If the interference can not be removed after taking the above measures, please use an EMG filter. The ECG waveform recorded at this time will be slightly attenuated.

7.4 Baseline drift

- Is the electrode installation stable?
- Is the connection of lead cables or electrodes reliable?
- Are the electrodes and patient's skin daubed with enough conductive paste?
- Is it caused by the patient's movement or breathing?
- Are the electrodes or leads in bad condition?

Warning: if the interference can not be removed after taking the above measures, please use a baseline filter

7.5 Troubleshooting list

Problem	Cause	Solution
Too large interference, disorderly waveform	<ol style="list-style-type: none"> 1. Grounding cable is not connected reliably. 2. Lead cables are not connected reliably. 3. There is AC interference. 4. Patient is nervous and is speaking. 	<ol style="list-style-type: none"> 1. Check the power cord and lead cables. 2. Let the patient prepare for the measurement.
Baseline burr	<ol style="list-style-type: none"> 1. AC interference is large. 2. Patient nervous and EMG interference is large. 	<ol style="list-style-type: none"> 1. Improve the environment. 2. If the bed is made of steel, replace it. 3. The power cable and lead cables are not parallel or too close to each other.
Irregular waveform, large up and down beeline figure	<ol style="list-style-type: none"> 1. Bad electrode conductivity. 2. Low battery. 3. Bad connection between electrodes and patient's skin. 4. Loose connection between lead cables and the device's plug. 5. Bad connection between electrodes and lead cables. 	<ol style="list-style-type: none"> 1. Use high quality alcohol. 2. Clean electrode and the skin under the electrode with alcohol. 3. Charge battery.
Baseline draft	<ol style="list-style-type: none"> 1. Low power 2. Patient is moving 	<ol style="list-style-type: none"> 1. Charge battery. 2. Keep patient still.
Unclear waveform	<ol style="list-style-type: none"> 1. Low battery. 2. The printer head surface is dirty. 3. Thermal paper problem. 	<ol style="list-style-type: none"> 1. Charge the battery. 2. Cut off the power, clean the printer head with alcohol, air dry. 3. Replace the thermal print paper with specified one.

8. MAINTENANCE

8.1 Battery

The device is designed with built-in fully sealed and maintenance-free rechargeable lithium battery, also equipped with perfect auto-charging and discharging monitor system. When the device is connected to AC power supply, the battery will be charged automatically. Battery status will be displayed on the right edge of the LCD screen in powering on state, as shown in the table. After it is completely discharged, the battery needs 5 hours to charge to 90% and 5.5 hours to charge to full capacity.

No.	Icon	Description
a		The battery status is unknown, generally appears 1 min after turning on.
b		Using AC power supply and the battery is full or no battery is in the device.
c		Using battery, battery is full
d		Using battery, battery is at 3/4 capacity
e		Using battery, battery is at 1/2 capacity
f		Using battery, battery is at 1/4 capacity
g		Using battery and the battery is low. It is recommended to charge it before use or switch to AC power supply.

Note: When charging the battery, the displayed status of the battery level switches between icon f to c.

The device can continuously print for 1.5 hours or work for more than 4 hours in standby mode when battery is completely charged. When the device is powered by battery, a battery icon will be displayed on the LCD screen, showing the battery capacity in 5 modes. When the battery capacity is too low for the device to operate, the device will turn off automatically to avoid permanent damage to the battery.

Note: The above is obtained by printing demo waveform under the test environment of temperature at 25°C, speed of 25 mm/s and gain of 10 mm/V. In actual use, the operation time may be shortened due to operation conditions and environment.

The battery should be recharged in time after discharged completely. If not used for a long period, the battery should be recharged every 3 months, which can extend the life of the battery.

When the battery can not be recharged or works no more than 10 minutes after fully charged, please replace the battery.

Warning:

- Do not try to disassemble the sealed battery without permission. The replacement of the battery shall be carried out by professional maintenance personnel authorized by our company and the same model of rechargeable battery provided by our company should be used.
- Do not touch the positive and negative terminals of the battery directly with wire, otherwise there is a danger of fire.
- Do not use the battery near fire sources or in environments where the temperature exceeds 60°C. Do not heat the battery, throw it into fire and keep away from water.
- Do not puncture, hammer or strike the battery or destroy it by other means, otherwise it will cause the battery to overheat, smoke, deform or burn.
- Keep away from the battery when it appears to be leaking or emitting an unpleasant smell. If the battery electrolyte leaks onto the skin or clothes, clean with water immediately. If the electrolyte accidentally enters your eyes, do not rub them, immediately clean with water and see a doctor.
- If the battery reaches its service life or the battery smells, is deformed, discolored or appears distorted, please stop using the battery and dispose of it according with local regulations.

8.2 Recording paper

In order to ensure the quality of the ECG waveform, please use the high-speed thermal recording paper supplied or specified by the company. If you use unspecified recording paper, the recorded ECG waveform may be blurred, faded and the paper feeding may not be smooth. This may even increase the wear of the device and shorten the service life of important parts such as the thermal print head. For information on how to purchase such recording paper, please contact your distributor or the company.

- When using recording paper, do not use recording paper with wax on the surface or in grey/black colors. Otherwise, the wax will stick to the heating part of the print head, resulting in abnormal work or damage the print head.

- High temperature, humidity and sunlight may cause the recording paper to change color. Please keep the recording paper in a dry and cool place.
- Please do not place the recording paper under fluorescent light for a long time, otherwise it will affect the recording effect.
- Please do not put the recording paper together with the PVC plastic, otherwise the color of the recording paper will change.
- Please use a recording paper with the specified dimensions. Recording paper that does not meet the requirements may damage the thermal print head or silicone rubber roller.

8.3 Maintenance after use

- Press  button to shutdown the device.
- Unplug the power cord and lead cabled. Hold the header of the plug to disconnect and do not pull the cable with excessive force.
- Clean the device and accessories, cover them to protect against dust.
- Store the device in a cool and dry place, avoid strong vibrations when moving.
- When cleaning the device, do not immerse it in the cleaning solution. Power supply must be cut off before cleaning. Use neutral detergents for cleaning. Do not use any detergent or disinfectant containing alcohol.

8.4 Lead cables and electrodes

The connectivity of the lead cable can be detected by the multimeter. Check whether each wire of the lead cable is in good contact according to the following table. The resistance of each wire from the electrode plug to the corresponding pin in the lead cable plug should be less than 10 Ω . The integrity of the lead cable must be checked regularly. Any lead wire damage will cause a false waveform of the corresponding lead or all leads on the ECG. The lead cable can be cleaned with neutral solvent. Do not use detergent or germicide containing alcohol (Please do not immerse the lead cables in liquid cleaning solution).

Note: The resistance of the lead cable with defibrillation-proof protection function is about 10 K Ω .

Mark	L	R	C1	C2	C3	C4	C5	C6	F	N
Pin position	10	9	12	1	2	3	4	5	11	14

Bending or knotting will shorten the service life of the lead cable. When using it, please straighten the lead cable first. The electrode should be well stored. After long time use, the surface of the electrode may oxidize and discolor due to corrosion and other factors, which may affect the signal acquisition. In this case, the electrode must be replaced.

8.5 Silicone rubber roller

The silicone rubber roller should be smooth and free of stains, otherwise it will affect the ECG recording effect. In order to remove the stains on the roller, please use a clean soft cloth dampened with a small amount of alcohol to wipe it along the longitudinal direction and scroll the roller in the paper conveying direction while wiping until it is clean.

8.6 Cleaning of thermal print head

Dirt and dust on the surface of the TPH can affect the clarity of the waveform. To clean the print head surface, open the paper compartment cover after turning off the device, use a clean and soft cloth dampened with alcohol to wipe the surface gently. For the residual stains on print head, moist it with a little alcohol first, then wipe with a soft cloth. Never use hard object to scratch the surface, otherwise the print head will be damaged. Wait until the alcohol has evaporated, then close the paper compartment cover. The print head should be cleaned at least once a month during normal use.

8.7 Disposal of product scrap

The disposal of packaging materials, waste battery and end-of-life device should obey the local laws and regulations and user should treat the scrapped products and materials properly according to the laws and regulations and try to support the classification and recycling work.

8.8 Others

Do not open the device enclosure to avoid electric shock danger.

The device associated circuit schematics and critical parts list are only available to authorized service station or maintenance personnel who is responsible for maintenance of the device.

The device belongs to measuring instrument. User should send the device to national designated inspection institution according to the requirements of the national metrological verification procedure. The device shall be inspected at least once per year and all accessories should be inspected and maintained regularly (at least once every six months).

9. PACKING LIST AND ACCESSORIES

9.1 Accompanying accessories

When the device is shipped from the factory, the intact packaging should contain the following contents:

Name	Quantity
Electrocardiograph	1 pc
Chest electrodes (suction cup/electrode slice)	1 set (6 pcs)
Limb electrodes (limb clip)	1 set (4 pcs)
ECG lead cable	1 pc
Potential equalization wire	1 pc
Power cord	1 pc
User manual	1 pc
Recording paper	1 pc

9.2 Notes

- Please follow the instructions on the package when opening the package.
- After unpacking, please check the accessories and accompanying documents in accordance with the packing list, then start inspecting the device.
- If the packaging content does not meet the requirement or the device does not work properly, please contact our company immediately.
- Please use the accessories provided by our company, otherwise the performance and safety of the device may be affected. If accessories provided by other company need to be used, please consult first with the after-sales service of the company.
- The packaging should be kept properly for future use in regular maintenance or device repair.

APPENDIX I ECG AUTOMATED MEASUREMENT AND INTERPRETATION GUIDE

1. Preface

The appendix describes the functions of the ECG automated measurement and automated interpretation. It explains the specific implementation methods, algorithm and formulas related to these two functions, as well as the content output by the automated measurement and automated interpretation.

According to the requirement of IEC60601-2-51:2003 *Medical electrical equipment - Part 2-51: Particular requirements for safety, including essential performance of recording and analysing single channel and multichannel electrocardiographs, Clause 50 Accuracy of operating data*, the appendix gives a description of verification process and results of the performance for automated measurement and automated interpretation.

2. Automated measurement parameters and Automated interpretation items

2.1 Measurement parameters

No.	Parameter	Unit
1	HR	bpm
2	PR-Interval	ms
3	P-duration	ms
4	QRS-duration	ms
5	T-duration	ms
6	QT/QTc	ms
7	P/QRS/T electric axis	deg
8	R(V5)/S(V1)	mV
9	R(V5)+S(V1)	mV

2.2 Interpretation items

No.	Item
1	No abnormal
2	Sinus mode Bradycardia
3	Sinus mode Tachycardia
4	Left atrium Hypertrophy
5	Right atrium Hypertrophy
6	Dual atrium Hypertrophy
7	QRS low voltage
8	Cardiac electric axis normal
9	Left axis deviation
10	Right axis deviation
11	Completeness Right Bundle branch block
12	Completeness Left Bundle branch block
13	Completeness Right Bundle branch block
14	Completeness Left Bundle branch block
15	V1 shows RSR' type
16	Left anterior fascicular block
17	Left posterior fascicular block
18	Left ventricular hypertrophy
19	Right ventricular hypertrophy
20	I atrioventricular block
21	Early anteroseptal MI
22	Possible acute forepart anteroseptal MI
23	Old anteroseptal MI
24	Early anterior MI
25	Possible acute anterior MI
26	Old anterior MI
27	Early extensive anterior MI
28	Possible acute anterior MI
29	Old extensive anterior MI
30	Early apical MI
31	Acute apical MI
32	Old apical MI
33	Early anterolateral MI
34	Possible acute anterolateral MI
35	Old anterolateral MI

No.	Item
36	Early high lateral MI
37	Possible acute high lateral MI
38	Old high lateral MI
39	Early inferior MI
40	Possible acute inferior MI
41	Old infolateral MI
42	Early inferolateral MI
43	Possible acute interolateral MI
44	Old inferolateral MI
45	ST depression, mild anteroseptal myocardial ischemia
46	ST depression, mild anterol myocardial ischemia
47	ST depression, mild extensive anterior myocardial ischemia
48	ST depression, mild apical myocardial ischemia
49	ST depression, mild anterolateral myocardial ischemia
50	ST depression, mild high lateral myocardial ischemia
51	ST depression, mild inferior myocardial ischemia
52	ST depression, mild inferolateral myocardial ischemia
53	ST depression, anteroseptal myocardial ischemia
54	ST depression, anterol myocardial ischemia
55	ST depression, extensive anterior myocardial ischemia
56	ST depression, mild apical myocardial ischemia
57	ST depression, anterolateral myocardial ischemia
58	ST depression, high lateral myocardial ischemia
59	ST depression, inferior myocardial ischemia
60	ST depression, inferolateral myocardial ischemia

2.2 Intended use

The intended use of the Automated Measurement and Interpretation function is shown as below:

Application and diagnosis	To detect abnormalities of the heart of the human body, examination items refer to above description
Population	Teenagers and adults, age range: 12-87
Application site	Hospitals
Accuracy	The accuracy of this function is reflected by the balance performance of sensitivity and specificity
Others	This function does not generate any alarm when using, so it should be operated by professionals or trained personnel.


3. Algorithm description

This section describes the algorithm, formulas and judgment conditions for interpretation items related to functions of ECG automated measurement and automated interpretation.

The 12-lead sync ECG waveform passes through the filter (AC, EMG, DFT (if it has and is open)) into the module of automated measurement and automated interpretation.

The module of automated measurement and automated interpretation mainly includes process of finding cardiac impulse location, finding the beginnin/end of each wave, amplitude calculation, parameters calculation and interpretations judgment based on known parameters.

The workflow is shown as below:

3.1 Find the cardiac impulse location

A. Data preprocessing, obtain the absolute value trend of slope for each lead; then superimpose each absolute value, obtain the superimposed graph of absolute value of slope.

B. Smoothing filter the superimposed graph on average of width 80ms, obtain the analytical data source DDD.

C. dFin the cardiac impulse location, give an initial threshold value: If the distance from the previous qrs-complex to the current location is less than 150 ms, then give up the location. Otherwise, take the 1/4 of threshold value as a reference, find the beginning of qrs-complex within 100 ms before the current location.

When the value is less than the threshold value, it may be the end of qrs-complex. Take the $1/4$ threshold value as a reference, find the end of qrs-complex. If the found qrs-complex is wide, this qrs-complex shall be excluded. Otherwise, save the found qrs-complex.
D. Locate: after you have found the qrs-complex, search the max value point between the beginning point and end point in the ECG original data, mark the point as cardiac impulse location.

E. Dynamically threshold adjustment: after finding the cardiac impulse location, use the value at the cardiac impulse location for the dynamically adaptive adjustment of the threshold value. Define the threshold value as $1/3$ of the average of the nearest three cardiac impulses.

F. After finding the cardiac impulse location, compute the RR-interval and accumulate it with the previous RR-intervals, then count the number of accumulated RR-intervals.

G. Continue searching until the end of the data and calculate the global average value for RR-intervals at the same time.

3.2 Find the beginning/end for each wave

The beginning/end of qrs complex has been approached in above cardiac impulse locating process but it is mainly in order to assist to find the cardiac impulse location; in addition, the location is searched based on the slope threshold value, which is imprecise. Here, according to the found cardiac impulse location, the beginning/end of qrs-complex will be sought accurately. Name the cardiac impulse location as the peak of R-wave.

• Read data

a. Read one data of qrs-complex: take the peak of R-wave as reference, locate directly to the original ECG file, read a piece of data containing the qrs-complex.

b. Preprocessing: superimpose the absolute value of slope for 12-lead signals.

c. Use the preprocessed data to carry on the searching of QRS-complex, P-wave and T-wave as the followings.

d. Read the next data of qrs-complex, repeat step 2 and step 4 until the analyzing of all qrs-complex are finished.

• Find QRS-complex

a. Calculate the threshold value of S-wave: search the minimal value within 200 ms after the peak of R-wave, take the value that equals to minimal value plus 0.4, as the threshold value for finding the end of S-wave.

b. Find the beginning of Q-wave: take 0.5 as the threshold value, search forwards starting from R-wave, a point that's less than the threshold value, within 0 ms-200 ms before the peak of R-wave, which is the beginning of Q-wave.

c. Find the end of S-wave: search backwards starting from R-wave, a point that's less than the threshold value of the end of S-wave, within 0 ms-200 ms after the peak of R-wave, which is the end of S-wave.

• Find P-wave

a. Peak of P-wave: search the max value within 30 ms-100 ms before the beginning of Q-wave, temporarily mark the point as the peak of P-wave.

b. Find the end of P-wave: search the minimal value between the peak of P-wave and the beginning of Q-wave, the minimal value plus the 0.05 is the threshold value, use the threshold value to find the end of P-wave.

c. Find the beginning of P-wave: search the minimal value within 150ms before the peak of P-wave, the minimal value plus 0.06 is the threshold value, use the threshold value to find the beginning of P-wave.

d. If the found P-wave is still narrow, it means that P-wave doesn't exist.

• Find T-wave

a. Peak of T-wave: search the max value within 30 ms- 300 ms after the end of QRS-complex, save it as the peak of T-wave.

b. Threshold value of the beginning of T-wave: search the minimal value within 0 ms-100 ms after the end of QRS-complex, the minimal value plus $1/10$ of the peak value of T-wave is the threshold for finding the beginning of T-wave.

c. Threshold value of the end of T-wave: search the minimal value within 20 ms after the peak of T-wave, the minimal value plus $1/10$ of the peak value of T-wave is the threshold for finding the end of T-wave.

d. Find the beginning of T-wave: in the range between the minimal value in step b and the peak of T-wave, find a point that is less than the threshold value of the beginning of T-wave, the point is the beginning of T-wave.

e. Find the end of T-wave: in the range between the minimal value in step c and the peak of T-wave, find a point that is less than the threshold value of the end of T-wave, the point is the end of T-wave.

• Explanation of the equipotential segment

In searching the QRS-complex, this algorithm adopts the analysis method of superposition of the slopes for all leads, therefore, the equipotential segments before and after the QRS-complex are partly included in the start and end point of the QRS-complex. It depends on the number of leads containing equipotential segments. If there are more leads containing equipotential segments, the slope value will be smaller after superposition, so it is difficult to meet the threshold condition and only a small part of the equipotential segments is counted to the start and end points of the QRS-complex. On the contrary, if there are less leads containing equipotential segments, a large part of the equipotential segments will be counted to the start and end points of the QRS-complex. Anyway, the equipotential segments before and after the QRS-complex are partly included in the QRS-complex duration.

3.3. Amplitude measurement

After finding the position of each wave, i.e. the start and end points of P wave, QRS complex and T wave, use the following method to measure P, Q, R, S, ST and T waves of each lead.

• P-wave

Calculate the average value of the data 20ms before the start point of P wave and use this average value as the baseline of P wave. Find the max value between the start point and end point of P wave, the difference between the max value and the baseline would be the amplitude of P wave.

- Q/R/S wave

Calculate the average value of the data 10-30 ms before the start point of QRS complex and use this average value as the baseline QRS complex. Search boundary points that exceed the baseline from the start point of Q wave to the end of S wave. Each adjacent two boundary points forms a sub-wave. Determine whether each sub-wave is a recognizable minimum wave (see the definition below). If it is a recognizable minimum wave, first identify its direction. If it is above the QRS baseline, it is R wave, if it is below the baseline it is Q wave or S wave. Find the extreme value of this wave and the difference between the extreme value and the baseline is the amplitude of Q/R/S wave.

- ST segment

Take above baseline of QRS complex as the ST baseline. Calculate the differences between the ST baseline and the points at 40 ms and 60 ms after the end point of QRS complex and calculate the average value of these two differences, the average value is the amplitude of ST segment.

- T-wave

Calculate the average value of the data 20-50 ms after the end point of T-wave and average this value with the QRS baseline in 2, then use the result as the baseline of T-wave. Find the max value between the start point and end point of T-wave, the difference between the max value and the baseline would be the amplitude of T-wave.

- Recognition of minimum wave

The minimum wave can be recognized by the algorithm according to the requirement of IEC60601-2-51:2003 Medical electrical equipment - Part 2-51: Particular requirements for safety, including essential performance of recording and analysing single channel and multichannel electrocardiographs, Annex GG, Clause GG.5 Definition of waveforms, measurement of minimum waves. The wave that meet the following conditions is the minimum wave that can be recognized by the algorithm.

- The signal part under consideration shows clearly two opposite slopes with at least one turning point in between;
- The signal part under consideration deviates at least $30\mu\text{V}$ from the reference level for a duration of at least 6 ms;
- The minimum observable duration of wave under consideration is 12 ms and amplitude $\geq 30\mu\text{V}$.

3.4. Calculation after intervals determination

The following parameters are determined according to the requirement of IEC60601-2-51:2003 Medical electrical equipment - Part 2-51: Particular requirements for safety, including essential performance of recording and analysing single channel and multichannel electrocardiographs, Annex GG Definitions and rules for the measurement of ELECRTOCARDIOGRAMS.

No.	Parameter	Calculation
1	HR	$60/RR$ (1)
2	PR-INTERVAL	Q_s (2)- P_s (3)
3	P-duration	P_e (4) - P_s (3)
4	QRS-duration	S_e (5) - Q_s (2)
5	T-duration	T_e (7) - T_s (6)
6	QT	T_e (7) - Q_s (2)
7	QTc	QT/\sqrt{RR} (1)
8	P/QRS/T electric axis	Electric axis formula: $\arctan(2.0 \times (S_{III} + S_I), S_I \times \sqrt{3}) / \pi$ (8) P electric axis: S_{III} : voltage sum from the beginning point to point of P-wave on lead III S_I : voltage sum from the beginning point to point of P-wave on lead I QRS electric axis: S_{III} : voltage sum from the beginning point to end point of QRS-complex on lead III S_I : voltage sum from the beginning point to end point of QRS-complex on lead I T electric axis: S_{III} : voltage sum from the beginning point to the end point of T-wave lead III S_I : voltage sum from the beginning point to the end point of T-wave lead I
9	R(V5)	Height (voltage value) of R-wave on lead V5
10	S(V1)	Height (voltage value) of S-wave on lead V1

Note:

- (1) RR: RR-interval
- (2) Qs: beginning of the Q-wave
- (3) Ps: beginning of the P-wave
- (4) Pe: end of the P-wave
- (5) Se: end of the S-wave
- (6) Ts: beginning of the T-wave
- (7) Te: end of the T-wave
- (8) PI: 3.1415926

3.5 Interpretation judgments based on parameters

No.	Item	Rule of interpretation
1	No abnormal	No abnormalities detected
2	Sinus mode Bradycardia	Sinus P-wave. PR-interval between 110 ms-210 ms, HR \leq */min, general* \approx 50
3	Sinus mode Tachycardia	Sinus P-wave. PR-interval between 110ms-210ms, HR \geq */MIN, GENERAL* \approx 100
4	Left atrium Hypertrophy	P-wave of leads I, II, aVL shall meet the conditions: width increase of P-wave \geq 110 ms, or P-wave displays in double-peak type, value of peak to peak
5	Right atrium Hypertrophy	For leads I, II, aVF, amplitude of P-wave \geq 0.25 mV or P-wave is sharp.
6	Dual atrium Hypertrophy	For leads I, II, aVF, amplitude of P-wave \geq 0.25 mV and P-wave duration \geq 110 ms
7	QRS low voltage	Voltage of I-aVF limb leads $<$ 0.5 mV and voltage of V1-V6 chest leads $<$ 0.8 mV
8	Cardiac electric axis normal	QRS-axis between 30 to 90 degree
9	Left axis deviation	QRS-axis between -90 to -30 degree
10	Right axis deviation	QRS-axis between 120 to 180 degree
11	Completeness Right Bundle branch block	QRS-duration $>$ 120 ms, R-wave of lead V1 or VR is wide (width of R wave $>$ 80 ms)
12	Completeness Left Bundle branch block	QRS-duration $>$ 120ms, R-wave of lead V5 or V6 is wide
13	No Completeness Right Bundle branch block	QRS-duration $<$ 120 ms, R-wave of lead V1 or VR is wide.
14	No Completeness Left Bundle branch block	QRS-duration $<$ 120 ms, R-wave of lead V15 or V6 is wide (width of R-wave $>$ 80 ms)
15	V1 shows RSR' type	QRS-complex of lead V1 is RSR' type
16	Left anterior fascicular block	QRS-duration $<$ 110 ms, QRS axis $<$ -30 degree, lead I and lead aVL are qR type, and Q-wave duration $<$ 20 ms, lead II, III and aVF are rS type.
17	Left posterior fascicular block	QRS-duration $<$ 110 ms, QRS-axis $>$ 90 degree, lead I and lead aVL are rS type, lead II, III and aVF are qR type and Q-wave of lead II and III $<$ 20 ms.
18	Left ventricular hypertrophy	R amplitude of lead I $>$ 1.5 mV, R amplitude of lead V5 $>$ 2.5 mV, R amplitude of lead aVL $>$ 1.2mV, R amplitude of lead aVF $>$ 2mV, R amplitude of lead V5 minus S amplitude of lead V1 $>$ 4 mV (male) or 3.5 mV (female).
19	Right ventricular hypertrophy	R amplitude of lead aVR $>$ 0.5mV, R amplitude of lead V1 $>$ 1mV, R amplitude of lead V1 minus S amplitude of lead V5 $>$ 1.2mV, R amplitude of lead V1 is larger than S amplitude, R amplitude of lead V5 is smaller than S amplitude, R amplitude of lead V5 is smaller than S amplitude.
20	I atrioventricular block	PQ interval $>$ 210 ms
21	Early anteroseptal MI	Early myocardial infarction change of leads V1, V2, V3 no change of leads V4, V5

22	Possible acute forepart anteroseptal MI	Acute myocardial infarction change of leads V1, V2, V3, no change of leads V4, V5
23	Old anteroseptal MI	Old myocardial infarction change of leads V1, V2, V3, no change of leads V4, V5
24	Early anterior MI	Early myocardial infarction
25	Possible acute anterior MI	Acute myocardial infarction change of leads V3, V4, V5 no change of leads V1, V2, V6.
26	Old anterior MI	Old myocardial infarction change of leads V3, V4, V5, no change of leads V1, V2, V6
27	Early extensive anterior MI	Early myocardial infarction change of leads V1, V2, V3, V4, V5
28	Possible acute extensive anterior MI	Acute myocardial infarction change of leads V1, V2, V3, V4, V5
29	Old extensive anterior MI	Old myocardial infarction change of leads V1, V2, V3, V4, V5
30	Early apical MI	Early myocardial infarction change of leads V4, V5, no change of leads V1, V2, V3
31	Acute apical MI	Acute myocardial infarction change of leads V4, V5, no change of leads V1, V2, V3
32	Old apical MI	Old myocardial infarction change of leads V4, V5 no change of leads V1, V2, V3
33	Early anterolateral MI	Early myocardial infarction change of leads I, aVL, V4, V5, V6
34	Possible acute anterolateral MI	Acute myocardial infarction change of leads I, aVL, V4, V5, V6
35	Old anterolateral MI	Old myocardial infarction change of leads I, aVL, V4, V5, V6
36	Early high lateral MI	Early myocardial infarction change of leads I, aVL, no change of leads II, III, aVF, V4, V5, V6
37	Possible acute high lateral MI	Acute myocardial infarction change of leads I, aVL, no change of leads II, III, aVF, V4, V5, V6
38	Old high lateral MI	Old myocardial infarction change of leads I, aVL, no change of leads II, III, aVF, V4, V5, V6
39	Early inferior MI	Early myocardial infarction change of leads II, III, aVF, no change of leads I, aVL
40	Possible acute inferior MI	Acute myocardial infarction change of leads II, III, aVF, no change of leads I, aVL
41	Old inferior MI	Old myocardial infarction change of leads II, III, aVF, no change of leads I, aVL
42	Early inferolateral MI	Early myocardial infarction change of leads I, II, III, aVL, aVF
43	Possible acute inferolateral MI	Acute myocardial infarction change of leads I, II, III, aVL, aVF
44	Old inferolateral MI	Old myocardial infarction change of leads I, II, III, aVL, aVF
45	ST depression, mild anteroseptal myocardial ischemia	Mild ST-segment depression of leads V1, V2, V3, and no change of leads V4, V5
46	ST depression, mild anterior myocardial ischemia	Mild ST-segment depression of leads V3, V4, V5, and no change of leads V1, V2, V6
47	ST depression, mild extensive anterior myocardial ischemia	Mild ST-segment depression of leads V1, V2, V3, V4, V5
48	ST depression, mild apical myocardial ischemia	Mild ST-segment depression of leads V4, V5 and no change of leads V1, V2, V3
49	ST depression, mild anterolateral myocardial ischemia	Mild ST-segment depression of leads I, aVL, V4, V5, V6
50	ST depression, mild high lateral myocardial ischemia	Mild ST-segment depression of leads I, aVL and no change of leads II, III, aVF, V4, V5, V6

51	ST depression, mild inferior myocardial ischemia	Mild ST-segment depression leads of II, III, aVF, and no change of leads I, aVL
52	ST depression, mild inferolateral myocardial ischemia	Mild ST-segment depression of leads I, II, III, aVL, aVF
53	ST depression, anteroseptal myocardial ischemia	Severe ST-segment depression of leads V1, V2, V3 and no change of leads V4, V5
54	ST depression, anterior myocardial ischemia	Severe ST-segment depression of leads V3, V4, V5 and no change of leads V1, V2, V6
55	ST depression, extensive anterior myocardial ischemia	Severe ST-segment depression of leads V1, V2, V3, V4, V5
56	ST depression, apical myocardial ischemia	Severe ST-segment depression of leads V4, V5 and no change of leads V1, V2, V3
57	ST depression, anterolateral myocardial ischemia	Severe ST-segment depression of leads I, aVL, V4, V5, V6
58	ST depression, high lateral myocardial ischemia	Severe ST-segment depression of leads I, aVL and no change of leads II, III, aVF, V4, V5, V6
59	ST depression, inferior myocardial ischemia	Severe ST-segment depression of leads II, III, aVF and no change of leads I, aVL
60	ST depression, inferolateral myocardial ischemia	Severe ST-segment depression of leads I, II, III, aVL, aVF

Note:

Early myocardial infarction: normal Q-wave, ST elevation or ST slope elevation.

Acute myocardial infarction: abnormal Q-wave, ST elevation or ST slope elevation.

Old myocardial infarction: abnormal Q-wave, no ST elevation.

Abnormal Q-wave:

For leads I, II, III, avR, avL, avF, V3, V4, V5, V6, voltage of Q-wave < -0.3 mV or 4 times of negative wave of Q-wave > voltage of R-wave and R'-wave and/or Q-duration > 40 ms.

For leads V1, V2, voltage of Q-wave < -0.08 mV and Q-duration > 10 ms.

ST elevation:

For leads I, II, III, avR, avL, avF, V4, V5, V6, the voltage of ST segment at 60 ms point > 0.1 mV and for leads V1, V2, V3, the voltage at 60 ms point > 0.3 mV.

ST slope elevation:

Voltage of ST segment at 20 ms point > = voltage of J point, voltage at 40 ms point > = the one at 20 ms, voltage at 60 ms point > = the one at 40 ms, with change of ST elevation.

4. Data sources and data preprocessing

4.1 Data sources

According to the requirement of IEC60601-2-51:2003 *Medical electrical equipment - Part 2-51: Particular requirements for safety, including essential performance, of recording and analysing single channel and multichannel electrocardiograph*, the CSE measurement database, CSE diagnostic database, CTS calibration database and customized data shall be used to evaluate the function of automated measurements and automated interpretations.

Verification	Database	Database items
Automated measurement	CTS database	CAL05000 CAL10000 CAL15000 CAL 20000 CAL20002 CAL20100 CAL20110 CAL20160 CAL20200 CAL202010 CAL20260 CAL20500 CAL30000 ANE20000 ANE20001 ANE20002
	CSE measurement database	MA_0001~MA0125
Automated interpretation	CSE diagnostic database	D_0001~D_1220
	Customized data	000001~000549

4.2 CTS introduction

The CTS computerized ECG conformance testing project was launched in 1989 by the European Union. This project laid the foundation for computerized ECG conformance testing service. Currently, about 20 types of waveform have been designed derived from the test signals having an infinite length, these signals are part of the CTS-ECG test database and have proven their effectiveness in a series of official tests. According to the requirement of *IEC60601-2-51:2003 Medical electrical equipment - Part 2-51: Particular requirements for safety, including essential performance of recording analysing single channel and multichannel electrocardiograph Clause 50.101.1*, 13 data (CAL05000 CAL10000 CAL15000 CAL 20000 CAL20002 CAL20100 CAL20110 CAL20160 CAL20200 CAL202010 CAL20260 CAL20500 CAL30000) are used in the automated parameters verification for this test.

4.3 CSE Introduction

The EU CSE (Common Standards for Quantitative Electrocardiography) ECG database contains 3-lead measurement database of collection1 and collection2, 12-lead measurement database of collection3 and collectio4 and a diagnostic database of collection5. In which, the 12-lead measurement database contains 250 groups of interference data; Diagnostic database contains 1220 cases of short-term ECG recording. The primary development purpose of using 12-lead or 15-lead is to evaluate the performance of the automatic ECG analyzer. In addition to the normal data, the database also includes clinically confirmed ECGs of variety cases, such as left ventricular hypertrophy, right ventricular hypertrophy, every part of myocardial infarction and ventricular hypertrophy accompanying myocardial infarction. The database has made a great contribution to the study of electrocardiology, which is, the CSE group published a report on the recommended standard for general ECG measurements based on the investigation and study of the database, which has been widely recognized by the world.

CSE database diagnostic items:

Item	Number
Normal	382
Left ventricular hypertrophy	183
Right ventricular hypertrophy	55
Biventricular hypertrophy	53
Anterior myocardial infarction	170
Inferior myocardial infarction	273
Complex myocardial infarction	104
Synthetical accuracy	1220

4.4 Customized data

4.4.1 Data description

Customized data	Description
Total recording number	549
Race	Asian
Converage of age, gender	Aged from 17 to 87, average age 57.23, standard deviation 21.32; 326 male, average age 55.54, standard deviation 19.81; 223 female, average age 59.70, standard deviation 22.63.
Sampling data	12-lead ECG data (I, II, III, AVR, AVL, AVF, V1, V2 V3, V4, V5, V6), sampling frequency of each channel: 1 kHz, amplitude quantization: 2.4 μ V/LSB
Remark	The interpretation conclusion of customized data is determined by the physician diagnostic results of cardiac catheterization and ultrasonic examination and the ECG judgment result in physical examination, the details as below: 1) Normal ECG: Determined by the diagnostic result the judged as normal in caridial catheterization and ultrasonic examination and the result that judged as normal in physical examination. 2) Atrium hypertrophy: Determined by the diagnostic results of ultrasonic examination. 3) Myocardial infarction and myocardial ischemia: Determined by the physician diagnostic results of cardiac catheterization. 4) Tachycardia, bradycardia, low voltage, axis: Determined by the diagnostic results of ultrasonic examination 5) Conduction block: Determined by the physician diagnostic results of cardiac catheterization. The standard of normal population in the customized database: physical examination is normal, no heart disease or other diseases that may affect cardiac functions of shape.

Note:
The heart abnormalities such as posterior myocardial ischemia, early posterior MI and old posterior MI are not included in the database. These abnormalities and other heart disorders not contained in above sheet won't be regarded as the judgment object for the verification of automated interpretation accuracy.

4.5. Data preprocessing

4.5.1 CTS preprocessing

The 16 cases (CAL05000, CAL10000, CAL15000, CAL20000, CAL20002, CAL20100, CAL20110, CAL20160, CAL20200, CAL20210, CAL20260, CAL20500, CAL30000, ANE20000, ANE20001, ANE20002) from CTS-ECG shall be processed for voltage conversion and frequency conversion for resampling as the applicable format in the system. Then cases will be imported to the device. After that, the verification of automated measurement parameters will be carried on.

4.5.2 CSE preprocessing

The cases (MA_0001-MA0125, D_0001~D_1220) from the CSE shall be processed for voltage conversion and frequency conversion for resampling as the applicable format in the system. Then cases will be imported to the device. After that, the case of MA_0001-MA0125 shall be used for the following verification of automated measurement parameters and the case of D_0001~D_1220 shall be used for the following verification of automated interpretation.


4.5.3 Customized data preprocessing

The customized initial case files shall be processed for voltage conversion and frequency conversion for resampling as the applicable format in the system. Then cases will be imported to the device. After that, the verification of automated interpretation will be carried on.

5. Process and Result of Verification

5.1 Verification of measurement function

• Verification and Process for CTS measurement database: (CAL05000, CAL10000, CAL15000, CAL20000, CAL20002, CAL20100, CAL20110, CAL20160, CAL20200, CAL20210, CAL20260, CAL20500, CAL30000, ANE20000, ANE20001, ANE20002) imported to the device shall be used to verify the automated measurement parameters.


• Verification and Process for CSE measurement database: Import the converted case files into the device, add appropriate database records, then waveform for all case files can be reviewed in the device, therefore the automated measurement parameters can be obtained.

Eliminate the cases existing obvious error for the diagnostic parameters (P-wave location is wrong) from the CSE database.

Make a comparison between the ECG analytical parameters (the beginning/end of P-wave, QRS-complex and T-wave) and the diagnostic parameters (the beginning/end of P-wave, QRS-complex and T-wave) provided by CSE database. Draw the two groups of waveform and mark the location of the beginning/end of P-wave, QRS-complex and T-wave corresponding to each case. The picture provides a visualized comparison, so the mean and standard deviation of the differences can be calculated. According to the requirement of IEC60601-2-21:2003 Medical electrical equipment - Part 2-51: Particular requirements for safety, including essential performance of recording and analysing single channel and multichannel electrocardiograph, the four largest deviations from the mean shall be eliminated before recalculation of mean and standard deviation of the differences.

Flow diagram of CSE measurement database verification process.

5.1.3 Verification results

- Accuracy of amplitude measurements

Calibration and analytical ECGs shall be used to measure the amplitude value, the summary as follows:

Amplitude	Mean difference (uV)	Standard deviation (uV)
P-wave	-1.70	5.72
Q-wave	7.51	18.07
R-wave	-18.05	21.70
S-wave	7.77	18.58
ST-segment	0.15	4.24
T-wave	-5.81	8.03

Note: In amplitude measurement, for large-amplitude ECG, such as CAL30000, it is necessary to adjust to 0.5 times the gain before testing.

- Accuracy of absolute interval and wave duration measurements

Calibration and analytical ECGs shall be used to measure the global interval and wave duration (including Q-wave, R-wave, S-wave), the summary as follows:

Interval and duration	Mean difference (ms)	Standard deviation (ms)
P-duration	-5.70	1.88
PQ-interval	-2.58	1.94
QRS-duration	-0.23	3.26
QT-interval	-6.70	4.37

- Accuracy of interval measurements on biological ECGs

CSE database shall be used to evaluate the accuracy of interval measurements on biological ECGs, the summary as follows:

Interval and duration	Mean difference (ms)	Standard deviation (ms)
P-duration	0.99	13.46
PR-interval	3.65	9.68
QRS-duration	-1.69	6.11
QT-interval	-2.32	20.69

- Stability of measurements against NOISE

The test is carrying on according to MA-series data (008, 011, 013, 014, 015, 021, 026, 027, 042, 061) in CSE database.

Global measurement parameters	Type of added NOISE	Disclosed differences	
		Mean (ms)	Standard deviation (ms)
P-duration	High frequency	-5.65	12.33
P-duration	Line frequency	-0.25	12.71
P-duration	Base-line	-4.90	33.15
QRS-duration	High frequency	-0.95	5.13
QRS-duration	Line frequency	1.35	4.71
QRS-duration	Base-line	-1.55	7.68
QT-interval	High frequency	-14.55	6.51
QT-interval	Line frequency	-8.55	20.73
QT-interval	Base-line	36.20	64.47

The biological ECGs are fed into the device in form of digital signals, then the measurement value can be obtained by calculation.

Test condition:


- a) without NOISE
- b) with 25 μ V high frequency
- c) with 50 μ V peak to valley 50Hz/60Hz sinusoidal line frequency NOISE
- d) with 1mV peak to valley 0.3Hz sinusoidal base-line NOISE

For each NOISE level above, the differences of measurements between the NOISE-free ECG and the ECGs with NOISE shall be determined. The two largest deviations from the mean shall be estimated before calculation of mean and standard deviation of differences.


5.2 Verification of interpretation function

5.2.1 Verification process

- CSE diagnostic database

• Customized database

5.2.2 Verification results

No.	Item	ECGs number	Sensitivity %	Specifity %	Positive predictive value %
1	No abnormal	585	92.01	79.16	97.38
2	Sinus mode Baradycardia	191	96.68	99.73	98.64
3	Sinus mode Tachycardia	78	97.44	96.49	96.90
4	Left atrium Hypertrophy	51	51.09	99.89	81.82
5	Right atrium Hypertrophy	43	42.64	99.66	50.00
6	Dual atrium Hypertrophy	22	93.58	99.14	60.19
7	QRS low voltage	5	96.37	99.36	63.25
8	Cardiac electric axis normal	733	98.36	89.13	98.79
9	Left axis deviation	168	98.65	89.40	98.18
10	Right axis deviation	107	98.23	88.99	94.90

11	Completeness Right Bundle branch block	28	97.00	89.50	95.45
12	Completeness Left Bundle branch block	32	97.73	89.65	91.43
13	No completeness Right Bundle branch block	41	96.86	89.83	82.35
14	No completeness Left Bundle branch block	47	94.68	89.93	89.66
15	V1 shows RSR' type	13	90.32	91.14	65.12
16	Left anterior fascicular block	26	91.43	93.25	71.11
17	Left posterior fascicular block	18	89.29	97.37	52.63
18	Left ventricular hypertrophy	236	41.37	92.65	70.36
19	Right ventricular hypertrophy	108	39.75	93.47	65.39
20	I atrioventricular block	13	94.58	91.67	80.64
21	Early anteroseptal MI	10	83.33	99.94	90.91
22	Possible acute forepart anteroseptal MI	27	16.67	98.73	91.89
23	Old anteroseptal MI	26	92.00	98.90	86.47
24	Early anterior MI	77	93.90	88.22	71.96
25	Possible acute anterior MI	10	80.00	99.72	44.44
26	Old anterior MI	13	24.00	99.66	50.00
27	Early extensive anterior MI	24	79.67	99.43	41.18
28	Possible acute extensive anterior MI	16	81.82	99.66	75.00
29	Old extensive anterior MI	30	90.91	88.05	37.04
30	Early apical MI	15	88.32	87.21	88.54
31	Acute apical MI	21	78.12	78.66	53.85
32	Old apical MI	19	79.63	89.94	80.00
33	Early anterolateral MI	36	77.51	79.96	83.33
34	Possible acute anterolateral MI	9	28.57	99.77	33.33
35	Old anterolateral MI	14	70.00	93.60	50.00
36	Early high lateral MI	16	79.65	95.78	80.42
37	Possible acute high lateral MI	8	81.60	99.94	85.71
38	Old high lateral MI	23	81.82	99.66	60.00
39	Early inferior MI	31	88.89	95.00	40.00
40	Possible acute inferolateral MI	29	11.11	99.94	50.00
41	Old inferior MI	101	96.07	99.24	93.44
42	Early inferolateral MI	73	98.77	96.82	75.94
43	Possible acute inferolateral MI	29	11.11	99.94	50.00
44	Old inferolateral MI	28	84.62	99.83	78.57
45	ST depression, mild anteroseptal myocardial ischemia	7	75.36	99.55	46.67
46	ST depression, mild anterior myocardial ischemia	5	81.24	99.94	33.33
47	ST depression, mild extensive anterior myocardial ischemia	13	79.83	99.13	53.59

48	ST depression, mild apical myocardial ischemia	17	76.97	99.14	43.13
49	ST depression, mild anterolateral myocardial ischemia	25	77.54	99.08	37.64
50	ST depression, mild high lateral myocardial ischemia	21	80.64	99.14	47.39
51	ST depression, mild inferior myocardial ischemia	12	79.73	99.60	55.16
52	ST depression, mild inferolateral myocardial ischemia	20	80.59	99.26	50.61
53	ST depression, anteroseptal myocardial ischemia	4	85.41	99.72	44.44
54	ST depression, anterior myocardial ischemia	12	87.66	98.58	34.85
55	ST depression, extensive anterior myocardial ischemia	7	84.78	98.04	67.75
56	ST depression, apical myocardial ischemia	18	79.95	99.14	55.12
57	ST depression, anterolateral myocardial ischemia	13	87.42	98.97	59.09
58	ST depression, high lateral myocardial ischemia	16	90.06	99.31	57.14
59	ST depression, inferior myocardial ischemia	12	89.88	99.13	40.08
60	ST depression, inferolateral myocardial ischemia	6	91.39	99.16	50.47

Sensitivity: probability that a "True sample" would be determined as certain "Item" by automated interpretation function;
 Specificity: probability that a "True unfit sample" would be determined as certain "Unfit item" by automated interpretation function;
 Positive predictive value: probability that a determine "Unfit item" is a "True unfit item".

APPENDIX II EMC Guidance and Manufacturer Declaration

Table 1:

Guidance and manufacturer's declaration - electromagnetic emission	
The Electrocardiograph is intended for use in the electromagnetic environment specified below. The purchaser or the user of the device should assure that it is used in such environment	
Emission test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class A
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable

Table 2:

Guidance and manufacturer's declaration - electromagnetic immunity		
The Electrocardiograph is intended for use in the electromagnetic environment specified below. The purchaser or the user of the device should assure that it is used in such environment		
Immunity test	IEC60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output line	±2 kV for power supply lines Not applicable

Surge IEC 61000-4-5	±1 kV lines to lines ±2 kV lines to earth	±1 kV lines to lines ±2 kV lines to earth
Voltage dips, short interruptions and voltage variations on power supply input lines	<5% UT (< 95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycle 70% UT (30% dip in UT) for 25 cycle <5% UT (>95% dip in UT) for 5 sec	<5% UT (< 95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycle 70% UT (30% dip in UT) for 25 cycle <5% UT (>95% dip in UT) for 5 sec
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m

Table 3:

Guidance and manufacturer's declaration - electromagnetic immunity		
The Electrocardiograph is intended for use in the electromagnetic environment specified below. The purchaser or the user of the device should assure that it is used in such environment		
Immunity test	IEC60601 test level	Compliance level
Conducted RF IEC61000-4-6	3V 0,15 MHz - 80 MHz 6V in ISM bands between 0,15 MHz - 80 MHz	3V 0,15 MHz - 80 MHz 6V in ISM bands between 0,15 MHz - 80 MHz
Radiated RF IEC61000-4-3	3 V/m 80 MHz - 2.7 GHz	3 V/m 80 MHz - 2.7 GHz
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.		
NOTE 2 These guidelines may not apply in all situations- Electromagnetic propagation is affected by absorption and reflection from structures, object and people.		
Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.		

Table 4:

Guidance and manufacturer's declaration - electromagnetic immunity							
The [Code SI] is intended for use in the electromagnetic environment specified below. The customer or the user of the [Code SI] should assure that it is used in such an environment							
	Test frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380-390	TETRA 400	Pulse modulation b) 18 Hz	1,8	0,3	27
	450	380-390	GMRS 460, FRS 460	FM c) ±5 kHz deviation 1 kHz sine	2	0,3	28
	710	704-787	LTE Band 13, 17	Pulse modulation b) 217 Hz	0,2	0,3	9
	745						
	780						
	810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	2	0,3	28
	870						
930							

Radiated RF IEC6100 0-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	2	0,3	28
	1845						
	1970						
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0,3	28
	5240	5100-5800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0,2	0,3	9
	5500						
5785							

NOTE: if necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3

- A) For some services, only the uplink frequencies are included.
- b) The carrier shall be modulated using a 50% duty cycle square wave signal.
- c) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:

$$E=6/d\sqrt{p}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m

WARNING

- Don't use near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Active medical devices are subject to special EMC precautions and they must be installed and used in accordance with guidelines.

Note:

- The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocation or re-orienting the equipment.
- When the device is disturbed, the data measured may fluctuate, please measure repeatedly or in another environment to ensure its accuracy.

	CONTEC MEDICAL SYSTEMS CO., LTD No. 112 Qinhuang West Street, Economic & Technical Development Zone, Qinhuangdao, Hebei Province, P.R.C.
	SHANGHAI INTERNATIONAL HOLDING CORP. GMBH (EUROPE) Eiffestrasse 80, 20537 Hamburg, Germany
	GRUPO R. QUERALTÓ S.A. Polígono Industrial El Pino, Calle Pino Albar, 24, 41016, Sevilla (Spain)